## **Product Datasheet**

# AbSelect (TM) Mouse Antibody Purification Kit 830-0005

Unit Size: 1 Test Store at 4C.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at: www.novusbio.com/830-0005

Updated 4/29/2024 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications
Submit a review at www.novusbio.com/reviews/destination/830-0005



### 830-0005

AbSelect (TM) Mouse Antibody Purification Kit

| Product Information |  |
|---------------------|--|
| Unit Size           | 1 Test   |
| Concentration       | Concentration is not relevant for this product. Please see the protocols for proper use of this product.   |
| Storage             | Store at 4C.   |
| Product Description |  |
| Description         | Commercially available antibodies often contain substances (e.g. BSA, glycine, tris and azide) that interfere in labeling reactions with enzymes, biotin, streptavidin or fluorophores. The AbSelect Mouse Antibody Purification System quickly removes these contaminants. It can also be used to purify mouse antibodies from crude samples such as ascites fluid.  The mouse antibody to be purified or cleaned up is in a volume of 0.1ml to 0.5ml. Up to 150ug of antibody can be purified in each run. The method involves |
|                     | capturing the mouse antibody on AbSelect Mouse resin. The resin has a high affinity for mouse IgG molecules. Once the mouse antibody has bound to the AbSelect Mouse resin, unwanted substances can be removed by simply washing the resin. The purified product is then eluted and neutralized. The AbSelect Mouse Antibody Purification System is not suitable for use with antibodies from other species. It should be noted that the binding strength of Bovine IgG to the AbSelect Mouse resin is negligible.               |
|                     | The components of the AbSelect kit are fully compatible with the Lightning-Link and Lightning-Link Rapid antibody labeling system (available separately), which allows the purified antibody to be rapidly conjugated with a wide variety of labels.   |
| Kit Components      | 1 or 3 vials of AbSelect Mouse resin, 1 vial of 10x Binding buffer, 1 vial of Washing buffer, 1 vial of Elution buffer, 1 vial of Neutralization buffer, 1 or 3 spin cartridge / collecting tube assemblies, 4 or 12 additional collecting tubes   |
| Notes               | AbSelect® is a registered trademark of Abcam and distributed by Novus Biologicals.  This product is for research use only and is not approved for use in humans or in clinical diagnosis. This product is guaranteed for 1 year from date of receipt and this statement overrides any mentioned guarantee period on the limitations section of this products datasheet. Please contact technical@novusbio.com with questions.  |





#### Novus Biologicals USA

10730 E. Briarwood Avenue Centennial, CO 80112 USA

Phone: 303.730.1950 Toll Free: 1.888.506.6887 Fax: 303.730.1966

nb-customerservice@bio-techne.com

#### **Bio-Techne Ltd**

19 Barton Lane Abingdon Science Park Abingdon, OX14 3NB, United Kingdom Phone: (44) (0) 1235 529449

Free Phone: 0800 37 34 15 Fax: (44) (0) 1235 533420 info.EMEA@bio-techne.com

#### **Bio-Techne Canada**

21 Canmotor Ave Toronto, ON M8Z 4E6 Canada

Phone: 905.827.6400 Toll Free: 855.668.8722 Fax: 905.827.6402

canada.inquires@bio-techne.com

#### **General Contact Information**

www.novusbio.com Technical Support: nb-technical@biotechne.com

Orders: nb-customerservice@bio-techne.com

General: novus@novusbio.com

#### Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/830-0005

Earn gift cards/discounts by submitting a publication using this product: www.novusbio.com/publications

