

# Product Datasheet

## 5X Antigen-Down Coating Buffer 1 NB900-81049

Unit Size: 100 ml

Store at 4C. Do not freeze.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NB900-81049](http://www.novusbio.com/NB900-81049)

Updated 2/5/2017 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NB900-81049](http://www.novusbio.com/reviews/destination/NB900-81049)



**NB900-81049****5X Antigen-Down Coating Buffer 1**

| <b>Product Information</b>         |                                                                                                                                                                                                                         |
|------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>                   | 100 ml                                                                                                                                                                                                                  |
| <b>Concentration</b>               | Please see the protocols for proper use of this product. If no protocol is available, contact technical services for assistance.                                                                                        |
| <b>Storage</b>                     | Store at 4C. Do not freeze.                                                                                                                                                                                             |
| <b>Product Description</b>         |                                                                                                                                                                                                                         |
| <b>Specificity/Sensitivity</b>     | 5X Antigen-Down Coating Buffer 1 (CB2)                                                                                                                                                                                  |
| <b>Product Application Details</b> |                                                                                                                                                                                                                         |
| <b>Applications</b>                | ELISA                                                                                                                                                                                                                   |
| <b>Recommended Dilutions</b>       | ELISA                                                                                                                                                                                                                   |
| <b>Application Notes</b>           | Add 1 part CB2 to 4 parts deionized water and mix for 5-15 minutes. Add the antigen to be coated to the 1X buffer and mix for 5-15 minutes. Coat the plate by pipetting the solution into the wells of the ELISA plate. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA

Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
[novus@novusbio.com](mailto:novus@novusbio.com)

### **Novus Biologicals Canada**

461 North Service Road West, Unit B37  
Oakville, ON L6M 2V5  
Canada

Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
[canada@novusbio.com](mailto:canada@novusbio.com)

### **Novus Biologicals Europe**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
[info@bio-techne.com](mailto:info@bio-techne.com)

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [technical@novusbio.com](mailto:technical@novusbio.com)  
Orders: [orders@novusbio.com](mailto:orders@novusbio.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Support products are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NB900-81049](http://www.novusbio.com/reviews/submit/NB900-81049)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

