# **Product Datasheet**

# SOD (Superoxide Dismutase) Activity Assay Kit (Colorimetric) KA0783

Unit Size: 1 Kit

Store at 4°C.

www.novusbio.com



technical@novusbio.com

**Publications: 4** 

Protocols, Publications, Related Products, Reviews, Research Tools and Images at: www.novusbio.com/KA0783

Updated 6/27/2023 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications Submit a review at www.novusbio.com/reviews/destination/KA0783



### **KA0783**

SOD (Superoxide Dismutase) Activity Assay Kit (Colorimetric)


<b>Product Information</b>	
Unit Size	1 Kit
Concentration	Concentration is not relevant for this product. Please see the protocols for proper use of this product.
Storage	Store at 4°C.

Product Description		
Description	Quality control test: Standard curve Superoxide Dismutase Activity Assay Kit is used for detecting SOD activity (% inhibition rate).	
Species	Rat, Mammal, Rabbit	
Kit Components	WST Solution, SOD Enzyme Solution, SOD Assay Buffer, SOD Dilution Buffer	
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.	
Assay Type	Colorimetric	
Suitable Sample Type	Blood Samples, Cells and Tissue	
Sample Volume	20 uL	

Product Application Details		
Applications	ELISA, Functional, Quantification	
Recommended Dilutions	ELISA, Functional, Quantification	

# **Images**

SOD (Superoxide Dismutase) Assay Kit (Colorimetric) [KA0783] - The standard curve for is for the purpose of illustration only and should not be used to calculate unknowns. A standard curve should be generated each time the assay is performed.



#### **Publications**

Thams S, Lowry E, Larraufie M et al. A stem cell-based screening platform identifies compounds that desensitize motor neurons to endoplasmic reticulum stress Am. J. Chin. Med. 2018-10-04 [PMID: 30446391] (WB, Human, Mouse)

Afroz R, Tanvir EM, Karim N et al. Sundarban Honey Confers Protection against Isoproterenol-Induced Myocardial Infarction in Wistar Rats. Biomed Res Int 2016-05-03 [PMID: 27294126]

Ibrahim M, Mikail MA, Ahmed IA et al. Comparison of the effects of three different Baccaurea angulata whole fruit juice doses on plasma, aorta and liver MDA levels, antioxidant enzymes and total antioxidant capacity. Eur J Nutr 2017-05-17 [PMID: 28516253]

Songstad NT, Kaspersen KH, Hafstad AD et al. Effects of High Intensity Interval Training on Pregnant Rats, and the Placenta, Heart and Liver of Their Fetuses. PLoS One. 2015-11-13 [PMID: 26566220]





# **Novus Biologicals USA**

10730 E. Briarwood Avenue Centennial, CO 80112 USA

Phone: 303.730.1950 Toll Free: 1.888.506.6887

Fax: 303.730.1966 novus@novusbio.com

# **Bio-Techne Ltd**

19 Barton Lane Abingdon Science Park Abingdon, OX14 3NB, United Kingdom Phone: (44) (0) 1235 529449

Free Phone: 0800 37 34 15 Fax: (44) (0) 1235 533420 info.EMEA@bio-techne.com

## **Bio-Techne Canada**

21 Canmotor Ave Toronto, ON M8Z 4E6 Canada

Phone: 905.827.6400 Toll Free: 855.668.8722 Fax: 905.827.6402

canada.inquires@bio-techne.com

# **General Contact Information**

www.novusbio.com

Technical Support: technical@novusbio.com

Orders: orders@novusbio.com General: novus@novusbio.com

#### Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/KA0783

Earn gift cards/discounts by submitting a publication using this product: www.novusbio.com/publications

