

Product Datasheet

Rabbit anti-Chicken IgG (H+L) Secondary Antibody [Rhodamine] NB120-6750

Unit Size: 1.5 mg

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB120-6750

Updated 5/14/2020 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB120-6750



NB120-6750

Rabbit anti-Chicken IgG (H+L) Secondary Antibody [Rhodamine]

| Product Information | |
|------------------------------------|---|
| Unit Size | 1.5 mg |
| Concentration | LYOPH mg/ml |
| Storage | Store at 4C in the dark. |
| Clonality | Polyclonal |
| Preservative | 0.01% Sodium Azide |
| Reconstitution Instructions | Reconstitute with 1.0 ml deionized water |
| Isotype | IgG |
| Conjugate | Rhodamine |
| Purity | Immunogen affinity purified |
| Buffer | Lyophilized from 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2, 10 mg/mL Bovine Serum Albumin (BSA) - Immunoglobulin and Protease free |

| Product Description | |
|--------------------------------|--|
| Host | Rabbit |
| Species | Chicken |
| Specificity/Sensitivity | Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Rabbit Serum and Chicken IgG. Some light chain cross reactivity may be observed against chicken immunoglobulins. |
| Immunogen | Chicken IgG whole molecule |

| Product Application Details | |
|------------------------------------|--|
| Applications | Western Blot, Fluorophore-linked immunosorbent assay, Flow Cytometry, Immunocytochemistry/Immunofluorescence |
| Recommended Dilutions | Western Blot, Flow Cytometry 1:500 - 1:2500, Immunocytochemistry/Immunofluorescence 1:1000 - 1:5000, Fluorophore-linked immunosorbent assay 1:10000 - 1:50000 |
| Application Notes | This product is designed for ICC/IF, FLISA and fluorescent western blotting. This product is also suitable for multiplex analysis, including multicolor imaging, utilizing various commercial platforms. |





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Secondary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB120-6750

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

