

Product Datasheet

Ethenocytidine Antibody (6F5) - BSA Free NB600-444

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB600-444

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB600-444



NB600-444

Ethenocytidine Antibody (6F5) - BSA Free

Product Information	
Unit Size	0.1 ml
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	6F5
Preservative	0.1% Sodium Azide
Isotype	IgG2 Lambda
Purity	Protein G purified
Buffer	PBS

Product Description	
Host	Mouse
Species	All Species
Specificity/Sensitivity	This is specific for ethenocytidine and ethenodecytidine.
Immunogen	Ethenocytidine

Product Application Details	
Applications	ELISA, Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	ELISA 1:100-1:2000, Immunocytochemistry/ Immunofluorescence 1:50-1:200
Application Notes	This Ethenocytidine antibody is useful for Immunocytochemistry/Immunofluorescence and ELISA.

Publications

Young, TL, Santella, RM. Development of techniques to monitor for exposure to vinyl chloride: monoclonal antibodies to ethenoadenosine and ethenocytidine. Carcinogenesis;9(4):589-92. 1988-04-01 [PMID: 3356066]



Procedures

Serum protocol for Ethenocytidine Antibody (NB600-444)

Ethenocytidine Antibody (6F5):

Competitive ELISA

I Coating of Plates

DNA coating: DNA is dissolved in PBS at appropriate concentration. 0.1 ml is added/well and plates put in 37 degrees Celsius incubator to evaporate overnight. Alternatively, plates can be coated with a 2-fold higher concentration of DNA for 2 hrs at 37 degrees Celsius then used. Column 1 is not coated. These well will not be used for the assay (no blocking, no antibody and no secondary antibody) but will have substrate added for blanking the reader. Plates are stored in the refrigerator.

Protein coating: Proteins are dissolved in PBS at the appropriate concentration. 0.1 ml is added/well and plates put in 37 degrees Celsius incubator to evaporate overnight. Column 1 is again not coated. Plates are stored in the refrigerator.

An alternate protein coating condition is to dissolve the protein in 0.1 M sodium carbonate buffer pH 9.6. 0.1 ml is added/well and the plates are refrigerated for several hours or overnight. They cannot be used after 3 days.

1 M solution 1.59 g Na₂CO₃ + 2.93 g NaHCO₃/100ml

II Assay

1. Label assay sheet and determine which rows are to be used. Row 1 (A-H) is not used; it will be used to blank the spectrophotometer. Avoid using the outer rows if possible (i.e. 12A-H, H 1-12 and A 1-12).
2. Wash plate with wash buffer containing PBS-Tween and NaN₃ 3 x on each side (right side up and upside down). Shake out onto paper towel.
3. Add 0.2 ml/well of 1% FCS in wash buffer to block non specific binding. Solution of FCS should be made fresh.
4. Incubate 1 hr.
5. Preparation of inhibitor series (during incubation of plate with FCS). Calculate appropriate concentrations to give desired fmol/well=fmol/0.05 ml. Make serial dilutions by adding PBS or CT DNA to tubes followed by competitor.
6. Prepare antibody in 1% FCS washing buffer.
7. At end of incubation period, shake out solution from plate and tap onto paper towel to dry.
8. Add 0.05 ml of competitor to each well followed by 0.05 ml of diluted antibody. Be sure to run all controls including zero (no competitor), minus Ab (no antigen specific antibody but secondary antisera) and positive and negative controls.
9. Incubate for 90 min at 37 degrees Celsius.
10. Wash the plate with washing buffer 3 times on each side. Tap onto paper towels.
11. Secondary antisera - Use goat anti-mouse IgG-alkaline phosphatase for monoclonals and anti rabbit for polyclonals. Dilute as appropriate and add 0.1 ml/well.
12. Incubate for 90 min at 37 degrees Celsius.
13. Wash with wash buffer 3 x each side. Tap onto paper towel.
14. Wash plate 2 times with 0.01 M diethanolamine using the was bottle and covering the well completely each time. Tap onto paper towel. This step removes phosphate buffer which inhibits alkaline phosphatase activity.
15. Prepare the substrate - 2 tablets 95 mg/tablet) Sigma 104 in 10 ml 1 M diethanolamine, pH 8.6. Final concentration 1 mg/ml. Avoid physical contact of skin with the tablets since skin contains alkaline phosphatase. Add

0.1 ml/well

16. Incubate at 37 degrees Celsius and read absorbance at 405 nm. The absorbance of the 0 fmol standard should be between 0.5 and 1. Values above 2 are not usable since the reader may not be linear in this range.

Rinse water - One liter of H₂O + 2 ml 10% NaN₃

Wash buffer - One liter of 1 x PBS + 500 ul Tween 20 + 2 ml 10% NaN₃

Blocking buffer - Wash buffer + 1% FCS





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB600-444

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NB600-444C	Ethenocytidine Antibody (6F5) [DyLight 650]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB600-444

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

