

Product Datasheet

Rabbit anti-Bovine IgA Heavy Chain Secondary Antibody NB770

Unit Size: 1 ml

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Publications: 4

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB770

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB770



NB770**Rabbit anti-Bovine IgA Heavy Chain Secondary Antibody**

Product Information	
Unit Size	1 ml
Concentration	1.0 mg/ml
Storage	Store at 4C. Do not freeze.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Purity	Immunogen affinity purified
Buffer	Phosphate Buffered Saline (PBS)

Product Description	
Host	Rabbit
Species	Bovine
Specificity/Sensitivity	By immunoelectrophoresis and ELISA this reacts specifically with bovine IgA. This may cross react with IgA from other species.
Immunogen	Bovine IgA

Product Application Details	
Applications	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 1:1000-1:10000, ELISA, Immunohistochemistry 1:200 - 1:2000, Immunocytochemistry/ Immunofluorescence 1:200 - 1:2000, Immunohistochemistry-Paraffin

Publications

Vidal RM, Montero DA, Del Canto F et al. Safety and Immunogenicity of a Chimeric Subunit Vaccine against Shiga Toxin-Producing Escherichia coli in Pregnant Cows International journal of molecular sciences 2023-02-01 [PMID: 36769094] (ELISA)

Villot C, Chen Y, Pedgerachny K et al. Early supplementation of Saccharomyces cerevisiae boulardii CNCM I-1079 in newborn dairy calves increases IgA production in the intestine at 1 week of age J. Dairy Sci. 2020-07-16 [PMID: 32684462] (IHC-P, Bovine)

Sardoy MC, Anderson DE, George A et al. Standardization of a method to detect bovine sperm-bound antisperm antibodies by flow cytometry. Theriogenology . [PMID: 22925638]

Details:
Citation using the FITC form of this antibody.

Bajer AA, Garcia-Tapia D, Jordan KR et al. Peripheral blood-derived bovine dendritic cells promote IgG1-restricted B cell responses in vitro. J Leukoc Biol 2003-01-01 [PMID: 12525567]





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-
techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB770

NBP3-20648-100ug	IgA Heavy Chain Antibody (IGHA/7567R)
------------------	---------------------------------------

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Secondary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB770

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

