

Product Datasheet

Goat F(ab')₂ anti-Rabbit IgG (H+L) Secondary Antibody [FITC] (Pre-adsorbed) NBP1-75338

Unit Size: 0.5 mg

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-75338

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-75338



NBP1-75338Goat F(ab')₂ anti-Rabbit IgG (H+L) Secondary Antibody [FITC] (Pre-adsorbed)

Product Information	
Unit Size	0.5 mg
Concentration	1.5 mg/ml
Storage	Store at 4C. Do not freeze.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide
Isotype	IgG
Conjugate	FITC
Purity	Affinity purified
Buffer	10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2, 1 % (w/v) BSA, Protease/IgG free
Product Description	
Description	Purity > 90% based on SDS-PAGE Small amounts of intact IgG may be present. . Goat serum was obtained from healthy animals of US origin and under the care of a registered veterinarian.
Host	Goat
Species	Rabbit
Reactivity Notes	Based on IEP, no reactivity is observed to non-immunoglobulin rabbit serum proteins, serum proteins from bovine, human, or mouse, and IgG from human or mouse
Specificity/Sensitivity	Based on IEP, this Goat F(ab') ₂ anti-Rabbit IgG (H+L) Secondary Antibody [FITC] (Pre-adsorbed) heavy gamma chains on rabbit IgG and light chains on all rabbit immunoglobulins. This antibody has been pre-adsorbed against bovine, human, or mouse IgG or serum proteins.
Immunogen	This Goat F(ab') ₂ anti-Rabbit IgG (H+L) Secondary Antibody [FITC] (Pre-adsorbed) was developed against purified rabbit IgG (H&L).
Product Application Details	
Applications	Flow Cytometry, Immunomicroscopy
Recommended Dilutions	Flow Cytometry 1:20-1:2000, Immunomicroscopy 1:20-1:2000
Application Notes	This conjugate is suitable for immunomicroscopy and flow cytometry.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Secondary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-75338

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

