## **Product Datasheet**

# Ovary Tissue Slides (Normal)- Paraffin NBP2-30190

Unit Size: 5 Slides

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

**Publications: 1** 

Protocols, Publications, Related Products, Reviews, Research Tools and Images at: www.novusbio.com/NBP2-30190

Updated 10/23/2024 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications Submit a review at www.novusbio.com/reviews/destination/NBP2-30190



## NBP2-30190

Ovary Tissue Slides (Normal)- Paraffin

<b>Product Information</b>		
Unit Size	5 Slides	
Concentration	Concentration is not relevant for this product. Please see the protocols for proper use of this product.	
Storage	Store at 4C. Do not freeze.	
Reconstitution Instructions	These slides are paraffin coated to prevent sample oxidization, it is recommended that slides are first de-paraffinized by baking at 60 degrees C for 1 hour in a vertical orientation prior to performing antigen retrieval procedures.	
Product Description		

<b>Product Description</b>	<b>Prod</b>	luct	Desc	ription
----------------------------	-------------	------	------	---------


Froduct Description		
Description	Layout: 1X1 Diameter: 4 Thickness:Single Tissue Slides Pos:A01 SI No:1 Age:34 Sex:F Organ:Ovary Tissue Status:Normal Species:Human	
Species	Human	
Lysate Type	Tissue	
Lysate Tissue Condition	Normal	

## **Product Application Details**

Applications	Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Immunohistochemistry, Immunohistochemistry-Paraffin

## **Images**

Hematoxylin & Eosin Stain: Ovary Tissue Slides (Normal) [NBP2-30190] - ovary, normal Age/Sex: 34/F



### **Publications**

Zsolt Szeles, Gábor L Petheő, Bence Szikora, Imre Kacskovics, Miklós Geiszt A novel monoclonal antibody reveals the enrichment of NADPH oxidase 5 in human splenic endothelial cells. Scientific reports 2023-11-01 [PMID: 37821487]





## **Novus Biologicals USA**

10730 E. Briarwood Avenue Centennial, CO 80112 USA

Phone: 303.730.1950 Toll Free: 1.888.506.6887

Fax: 303.730.1966

nb-customerservice@bio-techne.com

## **Bio-Techne Canada**

21 Canmotor Ave Toronto, ON M8Z 4E6 Canada

Phone: 905.827.6400 Toll Free: 855.668.8722 Fax: 905.827.6402

canada.inquires@bio-techne.com

## **Bio-Techne Ltd**

19 Barton Lane Abingdon Science Park Abingdon, OX14 3NB, United Kingdom Phone: (44) (0) 1235 529449

Free Phone: 0800 37 34 15 Fax: (44) (0) 1235 533420 info.EMEA@bio-techne.com

## **General Contact Information**

www.novusbio.com Technical Support: nb-technical@biotechne.com

Orders: nb-customerservice@bio-techne.com

General: novus@novusbio.com

### **Products Related to NBP2-30190**

NBP2-30290 Huma

Human Ovary Tissue MicroArray (Cancer)

#### Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Slides are guaranteed for 3 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-30190

Earn gift cards/discounts by submitting a publication using this product: www.novusbio.com/publications

