# **Product Datasheet**

# Kidney Tissue Slides (Chronic Pyelonephritis)- Paraffin NBP2-30268

Unit Size: 5 Slides

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at: www.novusbio.com/NBP2-30268

Updated 10/23/2024 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications Submit a review at www.novusbio.com/reviews/destination/NBP2-30268



#### NBP2-30268

Kidney Tissue Slides (Chronic Pyelonephritis)- Paraffin


Product Information	
Unit Size	5 Slides
Concentration	Concentration is not relevant for this product. Please see the protocols for proper use of this product.
Storage	Store at 4C. Do not freeze.
Reconstitution Instructions	These slides are paraffin coated to prevent sample oxidization, it is recommended that slides are first de-paraffinized by baking at 60 degrees C for 1 hour in a vertical orientation prior to performing antigen retrieval procedures.
Duadret Description	

	recommended that slides are first de-paraffinized by baking at 60 degrees C for 1 hour in a vertical orientation prior to performing antigen retrieval procedures.
<b>Product Description</b>	
Description	Layout: 1X1 Diameter: 4 Thickness:Single Tissue Slides Pos:A01 SI No:1 Age:49 Sex:M Organ:Kidney Pathology:Chronic Pyelonephritis Tissue Status:Diseased Species:Human
Species	Human
Lysate Type	Tissue
Lysate Tissue Condition	Chronic Pyelonephritis

Product Application Details	
Applications	Immunohistochemistry-Paraffin, Immunohistochemistry
Recommended Dilutions	Immunohistochemistry, Immunohistochemistry-Paraffin

## **Images**

Hematoxylin & Eosin Stain: Kidney Tissue Slides (Chronic Pyelonephritis) [NBP2-30268] - Age:49, Sex:M, Organ:Kidney, Pathology:Chronic Pyelonephritis, Tissue Status:Diseased





### **Novus Biologicals USA**

10730 E. Briarwood Avenue Centennial, CO 80112 USA

Phone: 303.730.1950 Toll Free: 1.888.506.6887

Fax: 303.730.1966

nb-customerservice@bio-techne.com

#### **Bio-Techne Canada**

21 Canmotor Ave Toronto, ON M8Z 4E6 Canada

Phone: 905.827.6400 Toll Free: 855.668.8722 Fax: 905.827.6402

canada.inquires@bio-techne.com

#### **Bio-Techne Ltd**

19 Barton Lane Abingdon Science Park Abingdon, OX14 3NB, United Kingdom Phone: (44) (0) 1235 529449

Free Phone: 0800 37 34 15 Fax: (44) (0) 1235 533420 info.EMEA@bio-techne.com

#### **General Contact Information**

www.novusbio.com Technical Support: nb-technical@bio-

techne.com

Orders: nb-customerservice@bio-techne.com

General: novus@novusbio.com

#### **Products Related to NBP2-30268**

NBP2-30220 Human Kidney Tissue MicroArray (Cancer)

#### Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Slides are guaranteed for 3 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-30268

Earn gift cards/discounts by submitting a publication using this product: www.novusbio.com/publications

