

Product Datasheet

PAP Pen Fine Tip NBP2-31111

Unit Size: 1 Pen

Store at room temperature.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-31111

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-31111



NBP2-31111

PAP Pen Fine Tip

| Product Information | |
|----------------------------|--|
| Unit Size | 1 Pen |
| Concentration | Please see the protocols for proper use of this product. If no protocol is available, contact technical services for assistance. |
| Storage | Store at room temperature. |

| Product Application Details | |
|------------------------------------|--|
| Applications | Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Frozen, In-situ Hybridization |
| Recommended Dilutions | Immunohistochemistry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Frozen, In-situ Hybridization |
| Application Notes | <p>Preparation of Working Solutions</p> <ol style="list-style-type: none"> 1. Ready to be used directly on wet slides. <p>Protocol Recommendations</p> <ol style="list-style-type: none"> 1. Hold the pen point up and press it with your fingertip to take blocked air out from the pen barrel to prevent overflow. 2. To start, press the tip of the pen straight down until fluid starts soaking the tip and then draw a circle around the specimen. |





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-
techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Support products are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-31111

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

