# **Product Datasheet**

# S100B Antibody (15F9NB) - BSA Free NBP2-45267

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

**Publications: 1** 

Protocols, Publications, Related Products, Reviews, Research Tools and Images at: www.novusbio.com/NBP2-45267

Updated 10/17/2023 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications Submit a review at www.novusbio.com/reviews/destination/NBP2-45267



### NBP2-45267

S100B Antibody (15F9NB) - BSA Free


, ,	
<b>Product Information</b>	
Unit Size	0.1 mg
Concentration	1.0 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	15F9NB
Preservative	0.02% Sodium Azide
Isotype	IgG2b Kappa
Purity	Protein G purified
Buffer	PBS

<b>Product Description</b>	
Host	Mouse
Gene ID	6285
Gene Symbol	S100B
Species	Human
Immunogen	Partial recombinant human S100B protein (amino acids 2-92) [Uniprot: P04271].


Product Application Details	
	Western Blot, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 2 ug/ml, Flow Cytometry 1:500, Immunohistochemistry 1:100 - 1:500, Immunocytochemistry/ Immunofluorescence 1:50-1:200, Immunohistochemistry-Paraffin 1:100 - 1:500

### **Images**


Western Blot: S100B Antibody (15F9NB) [NBP2-45267] - Total protein from human brain (1) and mouse brain (2) was separated on a 4-15% gel by SDS-PAGE, transferred to 0.2 um PVDF membrane for 30 minutes and blocked in 5% non-fat milk in TBST. The membrane was probed with 2 ug/ml anti-S100B in 1% milk and detected with an anti-mouse HRP secondary antibody using chemiluminescence.




Immunocytochemistry/Immunofluorescence: S100B Antibody (15F9NB) [NBP2-45267] - Analysis in A431 cells were fixed for 10 minutes using 10% formalin and then permeabilized for 5 minutes using 1X TBS + 0.5% Triton X-100. The cells were incubated with S100B (15F4NB) at a 1:100 dilution overnight at 4 degrees Celsius and detected with DyLight 488 (Green) at a 1:500 dilution. Actin was detected with Phalloidin 568 (Red) at a 1:200 dilution. Nuclei were detected with DAPI (Blue) at 2.0 ug/ml in 1X PBS. Cells were imaged using a 40X objective.




Immunohistochemistry-Paraffin: S100B Antibody (15F9NB) [NBP2-45267] - Analysis in human prostate tissue.




Flow Cytometry: S100B Antibody (15F9NB) [NBP2-45267] - An intracellular stain was performed on SK-MEL-28 cells with S100B [15F9NB] Antibody NBP2-45267PE (blue) and a matched isotype control (orange). Cells were fixed with 4% PFA and then permeablized with 0.1% saponin. Cells were incubated in an antibody dilution of 2.5 ug/mL for 30 minutes at room temperature. Both antibodies were conjugated to Phycoerythrin.




Western Blot: S100B Antibody (15F9NB) [NBP2-45267] - Total protein from human brain (1) and mouse brain (2) was separated on a 4-15% gel by SDS-PAGE, transferred to 0.2 um PVDF membrane for 30 minutes and blocked in 5% non-fat milk in TBST. The membrane was probed with 2 ug/ml anti-S100B in 1% milk and detected with an anti-mouse HRP secondary antibody using chemiluminescence.




Immunocytochemistry/Immunofluorescence: S100B Antibody (15F9NB) [NBP2-45267] - A431 cells were fixed for 10 minutes using 10% formalin and then permeabilized for 5 minutes using 1X TBS + 0.5% Triton X-100. The cells were incubated with S100B (15F9NB) at a 1:100 dilution overnight at 4 degrees Celsius and detected with DyLight 488 (Green) at a 1:500 dilution. Actin was detected with Phalloidin 568 (Red) at a 1:200 dilution. Nuclei were detected with DAPI (Blue) at 2.0 ug/ml in 1X PBS. Cells were imaged using a 40X objective.




Immunohistochemistry-Paraffin: S100B Antibody (15F9NB) [NBP2-45267] - Analysis of S100B (15F9NB) in human brain.




Immunohistochemistry-Paraffin: S100B Antibody (15F9NB) [NBP2-45267] - Analysis of S100B in human prostate.



Immunohistochemistry-Paraffin: S100B Antibody (15F9NB) [NBP2-45267] - Analysis in human brain tissue.



Flow Cytometry: S100B Antibody (15F9NB) [NBP2-45267] - Staining of A431 cells (1 x 10^6 cells/ml) with S100B antibody, clone 15F4.1D2 (Orange) stained at a dilution of 1:500. Detected with a GtxRb DyLight 488 secondary. Shown with the secondary control (Blue).



### **Publications**

Gao S, Tang Y, Sun W et al. 3D-bioprinted GelMA nerve guidance conduits promoted peripheral nerve regeneration by inducing trans-differentiation of MSCs into SCLCs via PIEZO1/YAP axis Materials Today Advances 2023-03-01 (ICC/IF, WB, Mouse, Rat)

#### **Procedures**

## Western Blot protocol for S100B Antibody (NBP2-45267)

- 1. Perform SDS-PAGE on samples to be analyzed, loading 25 ug of total protein per lane.
- 2. Transfer proteins to membrane according to the instructions provided by the manufacturer of the membrane and transfer apparatus.
- 3. Stain according to standard Ponceau S procedure (or similar product) to assess transfer success, and mark molecular weight standards where appropriate.
- 4. Rinse the blot.
- 5. Block the membrane using standard blocking buffer for at least 1 hour.
- 6. Wash the membrane in wash buffer three times for 10 minutes each.
- 7. Dilute anti-S100B primary antibody in blocking buffer and incubate 1 hour at room temperature.
- 8. Wash the membrane in wash buffer three times for 10 minutes each.
- 9. Apply the diluted HRP conjugated secondary antibody in blocking buffer (as per manufacturers instructions) and incubate 1 hour at room temperature.
- 10. Wash the blot in wash buffer three times for 10 minutes each (this step can be repeated as required to reduce background).
- 11. Apply the detection reagent of choice in accordance with the manufacturers instructions.

Note: Tween-20 can be added to the blocking or antibody dilution buffer at a final concentration of 0.05-0.2%.

# Immunocytochemistry/Immunofluorescence protocol for S100B Antibody (NBP2-45267)

S100B Antibody (15F9NB):

Immunocytochemistry Protocol

Culture cells to appropriate density in 35 mm culture dishes or 6-well plates.

- 1. Remove culture medium and add 10% formalin to the dish. Fix at room temperature for 30 minutes.
- 2. Remove the formalin and add ice cold methanol. Incubate for 5-10 minutes.
- 3. Remove methanol and add washing solution (i.e. PBS). Be sure to not let the specimen dry out. Wash three times for 10 minutes.
- 4. To block nonspecific antibody binding incubate in 10% normal goat serum from 1 hour to overnight at room temperature.
- 5. Add primary antibody at appropriate dilution and incubate at room temperature from 2 hours to overnight at room temperature.
- 6. Remove primary antibody and replace with washing solution. Wash three times for 10 minutes.
- 7. Add secondary antibody at appropriate dilution. Incubate for 1 hour at room temperature.
- 8. Remove antibody and replace with wash solution, then wash for 10 minutes. Add Hoechst 33258 to wash solution at 1:25,0000 and incubate for 10 minutes. Wash a third time for 10 minutes.
- 9. Cells can be viewed directly after washing. The plates can also be stored in PBS containing Azide covered in Parafilm (TM). Cells can also be cover-slipped using Fluoromount, with appropriate sealing.
- \*The above information is only intended as a guide. The researcher should determine what protocol best meets their needs. Please follow safe laboratory



#### Immunohistochemistry-Paraffin protocol for S100B Antibody (NBP2-45267)

- 1. Deparaffinize the tissue sections by immersing the slides in Xylene with two changes for 10 min each. Sections should not get dried at any stage from this point.
- 2. Rehydrate the tissue sections by immersing the slides in decreasing grades of ethanol as follows:
- a. Immerse in 100% ethanol with 2 changes for 5 minutes each
- b. Immerse in 95% ethanol with 2 changes for 5 minutes each
- c. Immerse in 90% ethanol for 5 minutes
- d. Immerse in 70% ethanol for 5 minutes
- e. Immerse in 50% ethanol for 5 minutes
- f. Immerse in distilled water for 5 minutes
- 3. Antigen Retrieval (Microwave Method):
- a. Immerse the slides in a microwave compatible tray containing 10 mM Sodium Citrate buffer (pH 6.0) with 0.05% Tween 20.
- b. Boil the slides and maintain the sub-boiling temperature for 5 minutes in the microwave. Thereafter, take out the tray very carefully and cool it at room temperature (RT) for about 30 minutes.
- c. Wash the slides 3 times, 3 minutes each by immersing them in TBST (Tris Buffered Saline having 0.05% Tween 20).
- Quenching of Endogenous Peroxidase:
- a. Incubate the slides in 3% hydrogen peroxide prepared in methanol for 15 minutes (at RT, in dark conditions).
- b. Wash the slides in TBST 3 times, 3 minutes each.
- Protein Blocking:
- a. Incubate the sections with background sniper solution at RT for 15 minutes (Biocare Medicals, USA).
- b. Wash the sections 3 times, 3 min each by immersing the slides in TBST.
- 6. Primary Antibody:
- a. Dilute the primary antibody at 5ug/ml concentration using PBS as a diluent.
- b. Incubate the sections with diluted primary antibody for 90 minutes at RT in a humidified chamber.
- c. Thereafter, wash the slides 4 times, 5 minutes each with TBST.
- 7. Probe (Secondary Reagent):
- a. Incubate with MACH 1 Mouse probe for 15 minutes at RT.
- b. Incubate for 30 min at room temperature with HRP-Polymer (Biocare Medical, USA).
- c. Wash the slides with TBST 4 times, 5 minutes each
- 8. Chromogen:
- a. Mix 32ul of DAB Chromogen with 1 ml of DAB substrate buffer (Biocare Medical, USA).
- b. Apply 200ul DAB mixture/section and incubate at RT in dark conditions (few seconds 5 minutes).
- c. As soon as an appropriate color develops, rinse the slides with deionized water (2-3 brief rinses).
- 9. Counter stain:
- a. Counter stain with Hematoxylin for 30 seconds (Vector Labs, USA).
- b. Wash in deionized water for 1-2 minutes to clear the extra stain.
- c. Incubate the slides in bluing solution or Scott's water twice for 2 minutes each time.
- 10. Dehydrate the sections in increasing grades of alcohols:
- a. 50% alcohol for 1 minute
- b. 70% for 1 minute
- c. 90% for 1 minute
- d. 95% for 1 minute
- e. 100% for 1 minute
- f. Xylene with 2 changes for 2 minutes each
- 11. Mount with DPX mount and cover-slip glass (Fisher Scientific, USA), carefully not allowing any air bubbles to enter.

NOTE:- This protocol is provided as a reference tool only. Depending upon the type of tissues /tissue processing and reagents employed, the end user will need to optimize the final conditions for achieving an expected staining.





## **Novus Biologicals USA**

10730 E. Briarwood Avenue Centennial, CO 80112 USA

Phone: 303.730.1950 Toll Free: 1.888.506.6887

Fax: 303.730.1966

nb-customerservice@bio-techne.com

# **Bio-Techne Ltd**

19 Barton Lane Abingdon Science Park Abingdon, OX14 3NB, United Kingdom

Phone: (44) (0) 1235 529449 Free Phone: 0800 37 34 15 Fax: (44) (0) 1235 533420 info.EMEA@bio-techne.com

#### **Bio-Techne Canada**

21 Canmotor Ave Toronto, ON M8Z 4E6 Canada

Phone: 905.827.6400 Toll Free: 855.668.8722 Fax: 905.827.6402

canada.inquires@bio-techne.com

### **General Contact Information**

www.novusbio.com Technical Support: nb-technical@biotechne.com

Orders: nb-customerservice@bio-techne.com

General: novus@novusbio.com

#### Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-45267

Earn gift cards/discounts by submitting a publication using this product: www.novusbio.com/publications

