Product Datasheet

8-Tube Magnetic Stand 265-1006

Unit Size: 1 Each

Store at room temperature.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at: www.novusbio.com/265-1006

Updated 2/6/2020 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications
Submit a review at www.novusbio.com/reviews/destination/265-1006



265-1006

8-Tube Magnetic Stand

0-1 ube Magnetic Stand	
Product Information	
Unit Size	1 Each
Concentration	Please see the protocols for proper use of this product. If no protocol is available, contact technical services for assistance.
Storage	Store at room temperature.
Product Description	
Description	The AbPure Magnetic Purification System is a convenient method to purify antibodies at a small scale with high recoveries. It quickly removes substances often contained in commercially available antibodies (e.g. BSA and other proteins, glycine, tris, and azide) that interfere with conjugation reactions. It can be easily adapted to your needs allowing you to purify from 20 to 200ug of antibody and to reach the desired final antibody concentration. The antibody to be purified or cleaned up should ideally be in a volume of 20ul to 1ml. The AbPure Magnetic Purification System has been designed to be compatible with our Lightning-Link, Thunder-Link PLUS, InnovaCoat GOLD and LATEX conjugation kits (available separately), which allow the purified antibody to be conjugated rapidly.
Notes	This product is manufactured by Expedeon Inc. and distributed by Novus Biologicals.
Product Application Details	

Notes	Biologicals.
Product Application Details	
Application Notes	This product complements the AbPure Magnetic Purification System (265-0200)





Novus Biologicals USA

10730 E. Briarwood Avenue Centennial, CO 80112 USA

Phone: 303.730.1950 Toll Free: 1.888.506.6887

Fax: 303.730.1966

nb-customerservice@bio-techne.com

Bio-Techne Ltd

19 Barton Lane Abingdon Science Park Abingdon, OX14 3NB, United Kingdom Phone: (44) (0) 1235 529449

Free Phone: 0800 37 34 15 Fax: (44) (0) 1235 533420 info.EMEA@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave Toronto, ON M8Z 4E6 Canada

Phone: 905.827.6400 Toll Free: 855.668.8722 Fax: 905.827.6402

canada.inquires@bio-techne.com

General Contact Information

www.novusbio.com Technical Support: nb-technical@biotechne.com

Orders: nb-customerservice@bio-techne.com

General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Support products are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/265-1006

Earn gift cards/discounts by submitting a publication using this product: www.novusbio.com/publications

