# **Product Datasheet**

# Golden Syrian Hamster IgG F(ab)2 Isotype Control NBP1-96782-1mg

Unit Size: 1 mg

Store at 4C short term. For extended storage, add an equal volume of glycerol, aliquot and store at -20C or below. Avoid repeated freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at: www.novusbio.com/NBP1-96782

Updated 8/29/2022 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications
Submit a review at www.novusbio.com/reviews/destination/NBP1-96782



# NBP1-96782-1mg

Golden Syrian Hamster IgG F(ab)2 Isotype Control

Golden Syrian Hamster IgG F(ab)2 Isotype Control	
Product Information	
Unit Size	1 mg
Concentration	1 mg/ml
Storage	Store at 4C short term. For extended storage, add an equal volume of glycerol, aliquot and store at -20C or below. Avoid repeated freeze-thaw cycles.
Preservative	0.01% Sodium Azide
Isotype	IgG
Purity	Ion exchange chromatography
Buffer	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Product Description	
Description	Native Golden Syrian Hamster IgG, Fab2.
Host	Golden Syrian Hamster
Species	Golden Syrian Hamster
Specificity/Sensitivity	Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Golden Syrian Hamster IgG, anti-Golden Syrian Hamster IgG F(ab)2 and anti-Golden Syrian Hamster Serum. No reaction was observed against anti-Golden Syrian Hamster IgG F(c) or anti-Pepsin.
Product Application Details	
Applications	Western Blot, ELISA, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 1:100 - 1:2000, ELISA 1:100 - 1:2000, Immunohistochemistry-

Paraffin

#### **Procedures**

### Antibody Competition protocol for IgG F(ab)2 Isotype Control (NBP1-96782)

Materials

Assay Buffer: 50 mM Tris, pH 8.0

Recombinant Human Pro TNF alpha Fusion Protein (rhPro-TNF-alpha ) (Catalog # 1012-PS)

Recombinant Human TACE/ADAM17 (rhTACE) (Catalog # 930-ADB)

Positive Control: Recombinant Human TNF-alpha (rhTNF-alpha) (Catalog # 210-TA)

Goat Anti-Human TNF alpha Biotinylated Antigen Affinity-purified Polyclonal Antibody (Catalog # BAF210)

SDS-PAGE followed by Western Blotting

- 1) Dilute rhPro-TNF-alpha to 0.2 mg/mL in Assay Buffer.
- 2) Dilute rhTACE to 0.2 mg/mL in Assay Buffer.
- 3) Prepare the following vials for a final volume of 20 uL.
  - a) 10 uL of rhPro-TNF-alpha at 0.2 mg/mL + 10 uL of rhTACE at 0.2 mg/mL
  - b) 10 uL of rhPro-TNF-alpha at 0.2 mg/mL + 10 uL of Assay Buffer (control-with incubation)
  - c) 10 uL of rhPro-TNF-alpha at 0.2 mg/mL + 10 uL of Assay Buffer (control-without incubation) Note: At this point the concentration of rhPro-TNF-alpha in the reaction tubes is 0.1 mg/mL.
- 4) Incubate vials at 37 degrees C overnight (except for the non-incubated control, store at less than or equal to 20 degrees C).
- 5) Stop the reactions by adding the reducing gel loading buffer for SDS-PAGE to all vials. Heat the samples at 100 degrees C for 3-5 minutes.
- a) Dilute the reaction tubes (and control without incubation) to 50 ng/15 uL in reducing sample buffer.
- 6) Prepare a sample of Positive Control at 5 ng/15 uL in reducing sample buffer. Heat the sample at 100 degrees C for 3-5 minutes.
- 7) Load the samples on a 15% gel.
  - a) 50 ng/lane (15 uL) of rhPro-TNF-alpha of the incubated reactions (including the control with incubation).
  - b) 50 ng/lane (15 uL) of rhPro-TNF-alpha of the control (without incubation).
  - c) 5 ng/lane (15 uL) of Positive Control (Catalog # 210-TA).
- 8) Follow SDS-PAGE/Western Blotting procedures.
- a) Use the Biotinylated Anti-human TNF-alpha /TNFSF1A antibody at 0.1 ug/mL. (Catalog # AF-210-NA may be used in place of # BAF210.)
- 9) Visually determine processing of rhPro-TNF-alpha to the mature form by rhTACE by comparing the incubated reactions to the Positive Control.<br/>
  | Caption:Support By Application | Caption:Support

https://www.novusbio.com/support/support-by-application.html

https://www.novusbio.com/spectraviewer





# **Novus Biologicals USA**

10730 E. Briarwood Avenue Centennial, CO 80112 USA

Phone: 303.730.1950 Toll Free: 1.888.506.6887

Fax: 303.730.1966

nb-customerservice@bio-techne.com

# **Bio-Techne Canada**

21 Canmotor Ave Toronto, ON M8Z 4E6 Canada

Phone: 905.827.6400 Toll Free: 855.668.8722 Fax: 905.827.6402

canada.inquires@bio-techne.com

## **Bio-Techne Ltd**

19 Barton Lane Abingdon Science Park Abingdon, OX14 3NB, United Kingdom Phone: (44) (0) 1235 529449

Free Phone: 0800 37 34 15 Fax: (44) (0) 1235 533420 info.EMEA@bio-techne.com

## **General Contact Information**

www.novusbio.com Technical Support: nb-technical@biotechne.com

Orders: nb-customerservice@bio-techne.com

General: novus@novusbio.com

#### Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Isotype Controls are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-96782

Earn gift cards/discounts by submitting a publication using this product: www.novusbio.com/publications

