

Product Datasheet

Flavivirus group antigen Antibody (D1-4G2-4-15 (4G2)) [Alexa Fluor® 647] NBP2-52709AF647

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Reviews: 1 Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-52709AF647

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-52709AF647



NBP2-52709AF647

Flavivirus group antigen Antibody (D1-4G2-4-15 (4G2)) [Alexa Fluor® 647]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	D1-4G2-4-15 (4G2)
Preservative	0.05% Sodium Azide
Isotype	IgG2a Kappa
Conjugate	Alexa Fluor 647
Purity	Protein A purified
Buffer	50 mM Sodium Borate

Product Description	
Host	Mouse
Species	Virus
Reactivity Notes	Dengue Virus, Zika Virus, West Nile Virus, Yellow Fever Virus, Flaviviridae.
Specificity/Sensitivity	This Flavivirus group antigen Antibody (D1-4G2-4-15 (4G2)) recognises flavivirus group specific antigens (Dengue virus, West Nile Virus, Japanese Encephalitis, Yellow Fever Virus, Zika virus etc). It binds to the fusion loop at the extremity of domain II of protein E.
Immunogen	This recombinant Flavivirus group antigen Antibody (D1-4G2-4-15 (4G2)) was prepared from Dengue Virus type 2 antigens.
Notes	Alexa Fluor (R) products are provided under an intellectual property license from Life Technologies Corporation. The purchase of this product conveys to the buyer the non-transferable right to use the purchased product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components, or any materials made using the product or its components, in any activity to generate revenue, which may include, but is not limited to use of the product or its components: (i) in manufacturing; (ii) to provide a service, information, or data in return for payment; (iii) for therapeutic, diagnostic or prophylactic purposes; or (iv) for resale, regardless of whether they are resold for use in research. For information on purchasing a license to this product for purposes other than as described above, contact Life Technologies Corporation, 5791 Van Allen Way, Carlsbad, CA 92008 USA or outlicensing@lifetech.com. This conjugate is made on demand. Actual recovery may vary from the stated volume of this product. The volume will be greater than or equal to the unit size stated on the datasheet.

Product Application Details	
Applications	Western Blot, ELISA, Flow Cytometry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin, CyTOF-ready, Immunofluorescence, Neutralization



Recommended Dilutions	Western Blot, Flow Cytometry, ELISA, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin, Neutralization, Immunofluorescence, CyTOF-ready
Application Notes	Optimal dilution of this antibody should be experimentally determined. Flavivirus group antigen antibody validated for IF from a verified customer review.

Images

IF assay in a high throughput format for screening Dengue virus inhibitors. HEK 293 cells in DMEM 5% FBS/1% PS at 5000 cells/well are infected in presence of inhibitors with Dengue 2 New Guinea (VR-1584) using 384-well collagen-coated plates. After 48 h incubation at 37C/5% CO₂, the supernatant is removed and the Alexa Fluor 647 conjugated Ab is added at 1:5250. Plates were read after O/N incubation. Columns 1-2: Cells control. 23-24, rows A-L: Cells + virus control. 23-24, rows M-P: Cells + virus control + control drug (Positive control). 3-12 and 13-22: Experimental compounds added in a dose response manner from high to low concentration (Columns 3-13 higher concentration, columns 12-22 lower concentration).



Publications

Ruofan Wang, Germán G. Gornalusse, Yeseul Kim, Urvashi Pandey, Florian Hladik, Lucia Vojtech Potent Restriction of Sexual Zika Virus Infection by the Lipid Fraction of Extracellular Vesicles in Semen *Frontiers in Microbiology* 2020-09-29 [PMID: 33133043]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-52709AF647

NBP1-96981AF647	Mouse IgG2a Kappa Isotype Control (M2AK) [Alexa Fluor® 647]
-----------------	---

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-52709AF647

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

