

# Product Datasheet

## Live-Dead Cell Staining Kit (Fluorometric) NBP2-54864

Unit Size: 96 Tests

Store at -20C.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-54864](http://www.novusbio.com/NBP2-54864)

Updated 4/29/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-54864](http://www.novusbio.com/reviews/destination/NBP2-54864)



**NBP2-54864****Live-Dead Cell Staining Kit (Fluorometric)**

<b>Product Information</b>	
<b>Unit Size</b>	96 Tests
<b>Concentration</b>	Concentration is not relevant for this product. Please see the protocols for proper use of this product.
<b>Storage</b>	Store at -20C.
<b>Product Description</b>	
<b>Description</b>	<p>Distinguish between live and dead cells within 20 min. Ready-to-Use Reagents. Detection method: Flow Cytometry (Ex/Em = 488/518 nm for stained live cells, Ex/Em = 488/615 nm for stained dead cells) Sample type: Adherent and suspension cells Applications: Stained live and dead cells can easily be visualized by fluorescence microscopy.</p> <p>This product is manufactured by Abcam and distributed by Novus Biologicals. (Abcam Catalog Number: ab65470)</p>
<b>Species</b>	Human, Mouse, Rat, Mammal
<b>Kit Components</b>	Solution A (1 mM Live-Dye), Solution B (1 mg/ml PI), Staining Buffer
<b>Suitable Sample Type</b>	Cell cultures (adherent and suspension cells)





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-54864](http://www.novusbio.com/reviews/submit/NBP2-54864)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

