

# Product Datasheet

## NorthernLights Streptavidin NL637 NL998-NOV

Unit Size: 0.5 ml

Store at 4C. Do not freeze.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NL998-NOV](http://www.novusbio.com/NL998-NOV)

Updated 9/11/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NL998-NOV](http://www.novusbio.com/reviews/destination/NL998-NOV)



**NL998-NOV****NorthernLights Streptavidin NL637**

<b>Product Information</b>	
<b>Unit Size</b>	0.5 ml
<b>Concentration</b>	Please see the protocols for proper use of this product. If no protocol is available, contact technical services for assistance.
<b>Storage</b>	Store at 4C. Do not freeze.
<b>Product Description</b>	
<b>Description</b>	Streptavidin NL637 is produced as the fluorochrome NL637 derivative of streptavidin. <b>Fluorochrome Absorption Maximum (nm) Emission Maximum (nm) NL637 637 658 Allophycocyanin 645 660 Alexa Fluor 647 650 668 Cy5 650 670</b> For use as a secondary developing reagent in immunofluorescence assays including immunohistochemistry, immunocytochemistry, flow cytometry, and other fluorescent immunoassays with biotinylated primary antibodies. <b>Precaution: Contains Sodium Azide as a preservative. Sodium Azide may react with lead and copper plumbing to form explosive metal azides. Flush with large volumes of water during disposal.</b>





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Support products are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NL998-NOV](http://www.novusbio.com/reviews/submit/NL998-NOV)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

