

Product Datasheet

FlexLISA (R) HRP Sandwich ELISA Optimization Kit 4200-0030

Unit Size: 1 Kit

Store at -20C.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/4200-0030

Updated 2/13/2020 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/4200-0030



4200-0030**FlexLISA (R) HRP Sandwich ELISA Optimization Kit**

Product Information	
Unit Size	1 Kit
Concentration	Concentration is not relevant for this product. Please see the protocols for proper use of this product.
Storage	Store at -20C.
Product Description	
Description	<p>The Expedeon Inc. FlexLISA kits are designed for the development and optimization of sandwich ELISA assays to detect the presence of any antigen with a high affinity for any antibody pair in complex samples (serum, plasma, urine etc.).</p> <p>Benefits of the FlexLISA Kit:</p> <ul style="list-style-type: none"> • Easy antibody labeling - uses world leading Lightning-Link technology • Uses 40 times less capture antibody - save on cost of materials • Use your antibodies of choice • Conjugate up to 3 capture and 3 detection antibodies - ideal for screening • Flexibility to run multiple conjugates on a plate simultaneously • Choice of HRP or AP for enzymatic detection • Rapid: 30 seconds hands-on time
Kit Components	3 x 10ug Lightning-Link® Horseradish Peroxidase (HRP), 3 x 10ug Lightning-Link® Streptavidin, 1 x black biotin pre-coated 96 stripwell plate (12 x 8 well strips) in a re-sealable foil pouch
Notes	<p>This product contains a black stripwell plate. For a clear stripwell plate, see 4200-0010.</p> <p>This product is manufactured by Expedeon Inc. and distributed by Novus Biologicals.</p>
Product Application Details	
Application Notes	A rapid one-step protocol is used in this kit, employing pre-blocked microtitre well strips coated with biotin. The FlexLISA assay is run by adding the sample and the antibody mix to the wells, incubating for 1 hour, washing and reading the assay plate. FlexLISA is ideal for quick and reliable ELISA assay development, antibody pair screening and ELISA assay optimization.



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom

Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/4200-0030

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

