

Product Datasheet

DAB Enhancing Solution H-2200-NB

Unit Size: 30 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H-2200-NB

Updated 1/26/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H-2200-NB



H-2200-NB**DAB Enhancing Solution**

Product Information	
Unit Size	30 ml
Concentration	Please see the protocols for proper use of this product. If no protocol is available, contact technical services for assistance.
Storage	Store at 4C in the dark.

Product Description	
Description	Features: DAB Enhancing Solution is a ready-to-use reagent for toning DAB stained sections. A ten second exposure to this solution darkens the reaction product. Only visible DAB precipitate is toned, so little amplification is observed with the enhancing reagent.
Notes	This product is manufactured by Vector Laboratories and distributed by Novus Biologicals.

Product Application Details	
Applications	Immunohistochemistry, Immunocytochemistry, Immunofluorescence
Recommended Dilutions	Immunohistochemistry, Immunofluorescence, Immunocytochemistry
Application Notes	Enhancement can take from 5 - 20 seconds and should be determined by the investigator. Usually a 10 - 15 second exposure to the Enhancing Solution is sufficient. The color may vary somewhat depending on the intensity of untuned DAB precipitate. Generally the toned DAB reaction product has a darker reddish brown color. Overexposed tissue will have increased background and may turn pink with time..

Images

DAB Enhancing Solution [H-2200-NB]





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Support products are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H-2200-NB

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

