

# Product Datasheet

## Colon adenocarcinoma tissue Tissue Slides (Adenocarcinoma)- Paraffin NBP2-61407

Unit Size: 5 Slides

Store at 4C under sterile conditions. Do not freeze.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-61407](http://www.novusbio.com/NBP2-61407)

Updated 10/23/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-61407](http://www.novusbio.com/reviews/destination/NBP2-61407)



**NBP2-61407**

Colon adenocarcinoma tissue Tissue Slides (Adenocarcinoma)- Paraffin


Product Information	
Unit Size	5 Slides
Concentration	Concentration is not relevant for this product. Please see the protocols for proper use of this product.
Storage	Store at 4C under sterile conditions. Do not freeze.
Reconstitution Instructions	These slides are paraffin coated to prevent sample oxidization, it is recommended that slides are first de-paraffinized by baking at 60 degrees C for 1 hour in a vertical orientation prior to performing antigen retrieval procedures.

Product Description	
Description	Layout: 1X1 Diameter: 4 Thickness:Single Tissue Slides Pos:A01 SI No:1 Age:46 Sex:F Organ:Colon Pathology:Adenocarcinoma Tissue Status:Metastasis to liver Species:Human
Species	Human
Lysate Type	Tissue
Lysate Tissue Condition	Adenocarcinoma
Lysate Protein State	Moderately Differentiated

Product Application Details	
Applications	Hematoxylin and Eosin Stain, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Immunohistochemistry, Immunohistochemistry-Paraffin, Hematoxylin and Eosin Stain

**Images**

Hematoxylin & Eosin Stain: Colon adenocarcinoma tissue Tissue Slides (Adenocarcinoma) [NBP2-61407]





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Slides are guaranteed for 3 months from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-61407](http://www.novusbio.com/reviews/submit/NBP2-61407)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

