

Product Datasheet

DAB/Ni Peroxidase Substrate Kit SK-4100-NB

Unit Size: 1 Kit

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/SK-4100-NB

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/SK-4100-NB




SK-4100-NB**DAB/Ni Peroxidase Substrate Kit**

Product Information	
Unit Size	1 Kit
Concentration	Concentration is not relevant for this product. Please see the protocols for proper use of this product.
Storage	Store at 4C in the dark.
Product Description	
Description	Features: Greater sensitivity than conventional substrates. Ideal for IHC/ICC/ISH and blots. Heat Stable Ideal for single and multiple labeling Stock solutions supplied in convenient dropper bottles promoting ease of handling. Includes nickel to produce a gray-black rather than brown reaction product.
Species	Human
Kit Components	Buffer Stock Solution (6 mL), DAB Stock Solution (6 mL), Nickel Solution (6 mL), Hydrogen Peroxide Solution (6 mL)
Notes	This product is manufactured by Vector Laboratories and distributed by Novus Biologicals.
Product Application Details	
Applications	Western Blot, Immunohistochemistry, Immunohistochemistry-Paraffin, In-situ Hybridization, Enzyme-linked Immunosorbent Spot, Immunocytochemistry
Recommended Dilutions	Western Blot, Immunohistochemistry, Immunohistochemistry-Paraffin, In-situ Hybridization, Enzyme-linked Immunosorbent Spot, Immunocytochemistry




Images

Immunohistochemistry: DAB/Ni Peroxidase Substrate Kit [SK-4100-NB] - Tumor tissue section showing specific cytoplasmic cell staining (brown, DAB) with hematoxylin QS (blue) counterstain.



Immunohistochemistry-Paraffin: DAB/Ni Peroxidase Substrate Kit [SK-4100-NB] - Paraffin-embedded human prostate tissue stained using rabbit anti-prostate specific antigen (PSA) detected with VECTASTAIN Elite ABC Reagent and Vector DAB substrate.



DAB/Ni Peroxidase Substrate Kit [SK-4100-NB]





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/SK-4100-NB

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

