## **Product Datasheet**

### VectaFluor(TM) Excel Amp. anti-Rabbit IgG DL488 Ab Kit DK-1488-NB

Unit Size: 1 Kit

Store at 4C. Do not freeze.

www.novusbio.com

technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at: www.novusbio.com/DK-1488-NB

Updated 6/25/2024 v.20.1

# Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications Submit a review at www.novusbio.com/reviews/destination/DK-1488-NB



#### DK-1488-NB


VectaFluor(TM) Excel Amp. anti-Rabbit IgG DL488 Ab Kit

Product Information	
Unit Size	1 Kit
Concentration	Concentration is not relevant for this product. Please see the protocols for proper use of this product.
Storage	Store at 4C. Do not freeze.
Conjugate	DyLight 488
Product Description	
Description	VectaFluor Excel Amplified DyLight® 488 Anti-Rabbit IgG offers a non-biotin amplification method for fluorescence applications. This system uses a ready-to- use (R.T.U.) Amplifier Antibody, followed by a ready-to-use (R.T.U.) VectaFluor DyLight dye-labeled detection system. <b>Features:</b> High sensitivity Low background VectaFluor DyLight 488 labeled anti-goat IgG antibody for detection Excitation: 493 nm Emission: 518 nm Color: Green
Kit Components	Normal Horse Serum (2.5%) for blocking - 15 mL, Amplifier Antibody (Goat anti- rabbit IgG) - 15 mL, VectFluor DyLight labeled Horse anti-Goat IgG - 15 mL
Notes	This product is manufactured by Vector Laboratories and distributed by Novus Biologicals.
Product Application Details	
Applications	Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	Immunocytochemistry/ Immunofluorescence



#### Images


Immunofluorescence: VectaFluor(TM) Excel Amp. anti-Rabbit IgG DL488 Ab Kit [DK-1488-NB] - Colon: Cytokeratin detected with VectaFluor Excel Amplified Kit, Anti- Rabbit IgG, DyLight 488 (green). Mounted in VECTASHIELD(R) HardSet(TM) Antifade Mounting Medium.



VectaFluor(TM) Excel Amp. anti-Rabbit IgG DL488 Ab Kit [DK-1488-NB]



VectaFluor(TM) Excel Amp. anti-Rabbit IgG DL488 Ab Kit [DK-1488-NB]







#### Novus Biologicals USA

10730 E. Briarwood Avenue Centennial, CO 80112 USA Phone: 303.730.1950 Toll Free: 1.888.506.6887 Fax: 303.730.1966 nb-customerservice@bio-techne.com

#### **Bio-Techne Canada**

21 Canmotor Ave Toronto, ON M8Z 4E6 Canada Phone: 905.827.6400 Toll Free: 855.668.8722 Fax: 905.827.6402 canada.inquires@bio-techne.com

#### **Bio-Techne Ltd**

19 Barton Lane Abingdon Science Park Abingdon, OX14 3NB, United Kingdom Phone: (44) (0) 1235 529449 Free Phone: 0800 37 34 15 Fax: (44) (0) 1235 533420 info.EMEA@bio-techne.com

#### **General Contact Information**

www.novusbio.com Technical Support: nb-technical@biotechne.com Orders: nb-customerservice@bio-techne.com General: novus@novusbio.com

#### Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/DK-1488-NB

Earn gift cards/discounts by submitting a publication using this product: www.novusbio.com/publications

