Product Datasheet

VectaFluor(TM) Excel Amp. anti-Mouse IgG DL488 Ab Kit DK-2488-NB

Unit Size: 1 Kit

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at: www.novusbio.com/DK-2488-NB

Updated 6/25/2024 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications
Submit a review at www.novusbio.com/reviews/destination/DK-2488-NB



DK-2488-NB

VectaFluor(TM) Excel Amp. anti-Mouse IgG DL488 Ab Kit


VectaFluor(TM) Excel Amp. anti-Mouse IgG DL488 Ab Kit	
Product Information	
Unit Size	1 Kit
Concentration	Concentration is not relevant for this product. Please see the protocols for proper use of this product.
Storage	Store at 4C. Do not freeze.
Conjugate	DyLight 488
Product Description	
Description	VectaFluor Excel Amplified Fluorescent Staining System offers a non-biotin amplification method for fluorescence applications. This system uses a ready-to-use (R.T.U.) Amplifier Antibody, followed by a ready-to-use (R.T.U.) VectaFluor DyLight dye-labeled detection system. Features: Ready-to-use High sensitivity Low background VectaFluor DyLight 488 labeled anti-goat IgG antibody for detection Excitation: 493 nm Emission: 518 nm Color: Green
Species	Mouse
Kit Components	Normal Horse Serum (2.5%) for blocking - 15 mL, Amplifier Antibody (Goat anti-Mouse IgG) - 15 mL, VectaFluor DyLight labeled Horse anti-Goat IgG - 15 mL
Notes	This product is manufactured by Vector Laboratories and distributed by Novus Biologicals.
Product Application Details	
Applications	Immunocytochemistry/ Immunofluorescence

Product Application Details	
Applications	Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	Immunocytochemistry/ Immunofluorescence



Images


Immunofluorescence: VectaFluor(TM) Excel Amp. anti-Mouse IgG DL488 Ab Kit [DK-2488-NB] - Tonsil: Multi-Cytokeratin detected with VectaFluor Excel Amplified Kit, Anti-Mouse IgG, DyLight 488 (green). Mounted in VECTASHIELD(R) HardSet(TM) Antifade Mounting Medium with DAPI (blue).



VectaFluor(TM) Excel Amp. anti-Mouse IgG DL488 Ab Kit [DK-2488-NB]



VectaFluor(TM) Excel Amp. anti-Mouse IgG DL488 Ab Kit [DK-2488-NB]





Novus Biologicals USA

10730 E. Briarwood Avenue Centennial, CO 80112 USA

Phone: 303.730.1950 Toll Free: 1.888.506.6887

Fax: 303.730.1966

nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave Toronto, ON M8Z 4E6 Canada

Phone: 905.827.6400 Toll Free: 855.668.8722 Fax: 905.827.6402

canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane Abingdon Science Park Abingdon, OX14 3NB, United Kingdom Phone: (44) (0) 1235 529449

Free Phone: 0800 37 34 15 Fax: (44) (0) 1235 533420 info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com Technical Support: nb-technical@biotechne.com

Orders: nb-customerservice@bio-techne.com

General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/DK-2488-NB

Earn gift cards/discounts by submitting a publication using this product: www.novusbio.com/publications

