

Product Datasheet

VectaFluor(TM) Excel Amp. anti-Mouse IgG DL594 Ab Kit DK-2594-NB

Unit Size: 1 Kit

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/DK-2594-NB

Updated 6/25/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/DK-2594-NB



DK-2594-NB

VectaFluor(TM) Excel Amp. anti-Mouse IgG DL594 Ab Kit

Product Information	
Unit Size	1 Kit
Concentration	Concentration is not relevant for this product. Please see the protocols for proper use of this product.
Storage	Store at 4C. Do not freeze.
Conjugate	DyLight 594

Product Description	
Description	VectaFluor Excel Amplified Fluorescent Staining System offers a non-biotin amplification method for fluorescence applications. This system uses a ready-to-use (R.T.U.) Amplifier Antibody, followed by a ready-to-use (R.T.U.) VectaFluor DyLight dye-labeled detection system. Features: Ready-to-use High sensitivity Low background VectaFluor DyLight 594 labeled anti-goat IgG antibody for detection Excitation: 592 nm Emission: 617 nm Color: Red
Species	Mouse
Kit Components	Normal Horse Serum (2.5%) for blocking - 15 mL, Amplifier Antibody (Goat anti-Mouse IgG) - 15 mL, VectaFluor DyLight labeled Horse anti-Goat IgG - 15 mL
Notes	This product is manufactured by Vector Laboratories and distributed by Novus Biologicals.

Product Application Details	
Applications	Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	Immunocytochemistry/ Immunofluorescence




Images


VectaFluor(TM) Excel Amp. anti-Mouse IgG DL594 Ab Kit [DK-2594-NB]



Immunofluorescence: VectaFluor(TM) Excel Amp. anti-Mouse IgG DL594 Ab Kit [DK-2594-NB] - Tonsil: Multi-Cytokeratin detected with VectaFluor Excel Amplified Kit, Anti-Mouse IgG, DyLight 594 (red). Mounted in VECTASHIELD(R) HardSet(TM) Antifade Mounting Medium with DAPI (blue).



VectaFluor(TM) Excel Amp. anti-Mouse IgG DL594 Ab Kit [DK-2594-NB]





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/DK-2594-NB

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

