

# Product Datasheet

## BoNT-D Light Chain Antibody (695004) [Janelia Fluor® 669] FAB6037JF669

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB6037JF669](http://www.novusbio.com/FAB6037JF669)

Updated 5/8/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB6037JF669](http://www.novusbio.com/reviews/destination/FAB6037JF669)



**FAB6037JF669**

BoNT-D Light Chain Antibody (695004) [Janelia Fluor® 669]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	695004
Preservative	0.05% Sodium Azide
Isotype	IgG1
Conjugate	Janelia Fluor 669
Purity	Protein A or G purified from hybridoma culture supernatant
Buffer	50mM Sodium Borate

Product Description	
Host	Mouse
Species	<i>C. botulinum</i>
Reactivity Notes	<i>C. botulinum</i>
Specificity/Sensitivity	Detects <i>C. botulinum</i> BoNT-D Light Chain in direct ELISAs and Western blots. In direct ELISAs and Western blots, less than 5% cross-reactivity with recombinant Light Chains of BoNT-A, -B, -C1, -E, -F, -G or recombinant Heavy Chains of BoNT-A, -C1, or -D is observed.
Immunogen	<i>E. coli</i> -derived recombinant <i>C. botulinum</i> BoNT-D Light Chain Thr2-Ser428 Accession # P19321
Notes	Sold under license from the Howard Hughes Medical Institute, Janelia Research Campus.

Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot
Application Notes	Optimal dilution of this antibody should be experimentally determined.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB6037JF669**

---

6037-ZN-010	BoNT-D Light Chain
-------------	--------------------

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB6037JF669](http://www.novusbio.com/reviews/submit/FAB6037JF669)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

