

# Product Datasheet

## Hematoxylin Counterstain

### H-3401-500-NB

Unit Size: 500 ml

Store at room temperature.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/H-3401-500-NB](http://www.novusbio.com/H-3401-500-NB)

Updated 2/14/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/H-3401-500-NB](http://www.novusbio.com/reviews/destination/H-3401-500-NB)




**H-3401-500-NB**

## Hematoxylin Counterstain


<b>Product Information</b>	
<b>Unit Size</b>	500 ml
<b>Concentration</b>	Please see the protocols for proper use of this product. If no protocol is available, contact technical services for assistance.
<b>Storage</b>	Store at room temperature.
<b>Product Description</b>	
<b>Description</b>	<p>Vector Hematoxylin is based on Gill's formulation, it produces blue stained nuclei with crisp detail.</p> <p>Alcohol and mercury free  Progressive stain formula  Provides excellent color contrast with most commonly used peroxidase and alkaline phosphatase substrates  Suitable for use with non-aqueous and aqueous mountants</p>
<b>Notes</b>	This product is manufactured by Vector Laboratories and distributed by Novus Biologicals.
<b>Product Application Details</b>	
<b>Applications</b>	Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, In-situ Hybridization
<b>Recommended Dilutions</b>	Immunohistochemistry, Immunocytochemistry/ Immunofluorescence, In-situ Hybridization

**Images**

Hematoxylin Counterstain [H-3401-500-NB] - Cytokeratin (m), ImmPRESS Anti-Mouse Ig Kit, ImmPACT DAB (brown) substrate. Hematoxylin (blue) counterstain.



Hematoxylin Counterstain [H-3401-500-NB] - Human breast carcinoma: Ki67 (m), VECTASTAIN Universal Quick Kit with DAB (brown) substrate. Hematoxylin (blue) counterstain.



Hematoxylin Counterstain [H-3401-500-NB]





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-  
techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Support products are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/H-3401-500-NB](http://www.novusbio.com/reviews/submit/H-3401-500-NB)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

