

Product Datasheet

Parainfluenza Virus type 2 Antibody [DyLight 680] NB120-8355FR-0.2ml

Unit Size: 0.2 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB120-8355FR

Updated 9/10/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB120-8355FR



NB120-8355FR-0.2ml

Parainfluenza Virus type 2 Antibody [DyLight 680]

| Product Information | |
|------------------------------------|--|
| Unit Size | 0.2 ml |
| Concentration | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. |
| Clonality | Monoclonal |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG2a |
| Conjugate | DyLight 680 |
| Purity | Protein G purified |
| Buffer | 50mM Sodium Borate |
| Product Description | |
| Host | Mouse |
| Species | Virus |
| Reactivity Notes | This antibody is reactive to Parainfluenza virus 2, but not type 1 or type 3. |
| Specificity/Sensitivity | This reacts with F protein of parainfluenza virus type 2. Purity > 90 % (SDS-PAGE) |
| Immunogen | MA 104 cells, inoculated with parainfluenza 2 virus, strain II ALTB cc 2056. Virus was grown on monolayer of MA cells. When the cultures showed pronounced c.p.e. virus containing medium was purified from cell fragments by low speed centrifugation (5000 rpm, 40 min). Virus was pelleted by high speed centrifugation at 18000 rpm for 2.5 hours in the SW27 rotor for 20 min at + 4 oC, then the supernatant was collected and protein concentration was determined by Lowry method. |
| Notes | Dylight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries. This conjugate is made on demand. Actual recovery may vary from the stated volume of this product. The volume will be greater than or equal to the unit size stated on the datasheet. |
| Product Application Details | |
| Applications | ELISA, Immunocytochemistry/Immunofluorescence |
| Recommended Dilutions | ELISA, Immunocytochemistry/Immunofluorescence |



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-
techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB120-8355FR

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

