Product Datasheet

Oval Cell Marker Antibody (OC2-3C5) [mFluor Violet 500 SE] NBP1-18968MFV500

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com

technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at: www.novusbio.com/NBP1-18968MFV500

Updated 3/28/2024 v.20.1

Earn rewards for product reviews and publications.



Submit a publication at www.novusbio.com/publications Submit a review at www.novusbio.com/reviews/destination/NBP1-18968MFV500

NBP1-18968MFV500

Oval Cell Marker Antibody (OC2-3C5) [mFluor Violet 500 SE]

Product Information

Unit Size	0.1 ml
Concentration	This product is unpurified. The exact concentration of antibody is not quantifiable.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	OC2-3C5
Preservative	0.05% Sodium Azide
Isotype	lgG2a
Conjugate	mFluor Violet 500 SE
Purity	Tissue culture supernatant
Buffer	50mM Sodium Borate

Product Description

•	
Host	Rat
Species	Mouse
Marker	Hepatic periductal cell Marker
Specificity/Sensitivity	Specific for both proliferating and normal hepatic periductal and ductal cells.
Immunogen	Nonparenchymal cells from DDC-treated mice to induce the oval cell response.
Notes	mFluor(TM) is a trademark of AAT Bioquest, Inc. This conjugate is made on demand. Actual recovery may vary from the stated volume of this product. The volume will be greater than or equal to the unit size stated on the datasheet.

Product Application Details	
Applications	Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Frozen
Recommended Dilutions	Flow Cytometry, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Frozen
Application Notes	Optimal dilution of this antibody should be experimentally determined.

Images

Oval Cell Marker Antibody (OC2-3C5) [mFluor Violet 500 SE] - Vial of mFluor Violet 500 conjugated antibody. mFluor Violet 500 is optimally excited at 410 nm by the Violet laser (405 nm) and has an emission maximum of 501 nm.

	mFluor™∖	/iolet 500)
	111110101	10101000	
	LASER (nm)	FILTER	7
mFluor" Violet 500	LASER (nm) Violet (405)	FILTER 515/20	
th Fluor ¹ Violet 50 ⁰	and the second se]





Novus Biologicals USA

10730 E. Briarwood Avenue Centennial, CO 80112 USA Phone: 303.730.1950 Toll Free: 1.888.506.6887 Fax: 303.730.1966 nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave Toronto, ON M8Z 4E6 Canada Phone: 905.827.6400 Toll Free: 855.668.8722 Fax: 905.827.6402 canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane Abingdon Science Park Abingdon, OX14 3NB, United Kingdom Phone: (44) (0) 1235 529449 Free Phone: 0800 37 34 15 Fax: (44) (0) 1235 533420 info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com Technical Support: nb-technical@biotechne.com Orders: nb-customerservice@bio-techne.com General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-18968MFV500

Earn gift cards/discounts by submitting a publication using this product: www.novusbio.com/publications

www.novusbio.com

