# **Product Datasheet**

# Lambda Light Chain Antibody (HP6054) [mFluor Violet 450 SE] NBP2-33309MFV450

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at: www.novusbio.com/NBP2-33309MFV450

Updated 10/26/2023 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications
Submit a review at www.novusbio.com/reviews/destination/NBP2-33309MFV450



# NBP2-33309MFV450

Lambda Light Chain Antibody (HP6054) [mFluor Violet 450 SE]

| Lambda Light Chain Antibody (HP6054) [mFluor Violet 450 SE]   |  |
|---|--|
| Product Information   |  |
| 0.1 ml  |  |
| Please see the vial label for concentration. If unlisted please contact technical services.   |  |
| Store at 4C in the dark.  |  |
| Monoclonal  |  |
| HP6054  |  |
| 0.05% Sodium Azide  |  |
| IgG2a Kappa   |  |
| mFluor Violet 450 SE  |  |
| Protein A or G purified   |  |
| 50mM Sodium Borate  |  |
| Product Description   |  |
| Mouse   |  |
| 3537  |  |
| IGLC1   |  |
| Human   |  |
| B-Cell Marker   |  |
| This monoclonal antibody is specific to lambda light chain of immunoglobulin and shows no cross-reaction with lambda light chain or any of the five heavy chains. In mammals, the two light chains in an antibody are always identical, with only one type of light chain, kappa or lambda. The ratio of Kappa to Lambda is 70:30. However, with the occurrence of multiple myeloma or other B-cell malignancies this ratio is disturbed. Antibody to the lambda light chain is reportedly useful in the identification of leukemias, plasmacytomas, and certain non-Hodgkins lymphomas. Demonstration of clonality in lymphoid infiltrates indicates that the infiltrate is malignant. |  |
| Purified human Lambda Light Chain myeloma proteins coupled to polyaminostyren microbeads  |  |
| mFluor(TM) is a trademark of AAT Bioquest, Inc. This conjugate is made on demand. Actual recovery may vary from the stated volume of this product. The volume will be greater than or equal to the unit size stated on the datasheet.   |  |
| Product Application Details   |  |
| Western Blot, Flow Cytometry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin, CyTOF-ready  |  |
| Western Blot, Flow Cytometry, Immunohistochemistry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin, CyTOF-ready  |  |
| Optimal dilution of this antibody should be experimentally determined.  |  |
|   |  |



# **Images**

Lambda Light Chain Antibody (HP6054) [mFluor Violet 450 SE] [NBP2-33309MFV450] - Vial of mFluor Violet 450 conjugated antibody. mFluor Violet 450 is optimally excited at 406 nm by the Violet laser (405 nm) and has an emission maximum of 445 nm.





# **Novus Biologicals USA**

10730 E. Briarwood Avenue Centennial, CO 80112 USA

Phone: 303.730.1950 Toll Free: 1.888.506.6887

Fax: 303.730.1966 novus@novusbio.com

#### **Bio-Techne Ltd**

19 Barton Lane Abingdon Science Park Abingdon, OX14 3NB, United Kingdom Phone: (44) (0) 1235 529449

Free Phone: 0800 37 34 15 Fax: (44) (0) 1235 533420 info.EMEA@bio-techne.com

#### **Bio-Techne Canada**

21 Canmotor Ave Toronto, ON M8Z 4E6 Canada

Phone: 905.827.6400 Toll Free: 855.668.8722 Fax: 905.827.6402

canada.inquires@bio-techne.com

### **General Contact Information**

www.novusbio.com

Technical Support: technical@novusbio.com

Orders: orders@novusbio.com General: novus@novusbio.com

#### Products Related to NBP2-33309MFV450

NB7552

Goat anti-Mouse Lambda Light Chain Secondary Antibody [HRP]

#### Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-33309MFV450

Earn gift cards/discounts by submitting a publication using this product: www.novusbio.com/publications

