

# Product Datasheet

## Kappa Light Chain Antibody (SPM508) [Alexa Fluor™ Plus 594]

**NBP2-34396AFP594**

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-34396AFP594](http://www.novusbio.com/NBP2-34396AFP594)

Updated 4/12/2026 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-34396AFP594](http://www.novusbio.com/reviews/destination/NBP2-34396AFP594)



**NBP2-34396AFP594**

Kappa Light Chain Antibody (SPM508) [Alexa Fluor™ Plus 594]

| <b>Product Information</b>     |  |
|--------------------------------|--|
| <b>Unit Size</b>               | 0.1 ml   |
| <b>Concentration</b>           | Please see the vial label for concentration. If unlisted please contact technical services.  |
| <b>Storage</b>                 | Store at 4C in the dark.   |
| <b>Clonality</b>               | Monoclonal   |
| <b>Clone</b>                   | SPM508   |
| <b>Preservative</b>            | 0.05% Sodium Azide   |
| <b>Isotype</b>                 | IgG1 Kappa   |
| <b>Conjugate</b>               | Alexa Fluor Plus 594   |
| <b>Purity</b>                  | Protein A or G purified  |
| <b>Buffer</b>                  | 50mM Sodium Borate   |
| <b>Product Description</b>     |  |
| <b>Host</b>                    | Mouse  |
| <b>Gene ID</b>                 | 3514   |
| <b>Gene Symbol</b>             | IGKC   |
| <b>Species</b>                 | Human, Rat (Negative)  |
| <b>Reactivity Notes</b>        | Does not react with Rat.   |
| <b>Marker</b>                  | B-Cell Marker  |
| <b>Specificity/Sensitivity</b> | This monoclonal antibody is specific to kappa light chain of immunoglobulin and shows no cross-reaction with lambda light chain or any of the five heavy chains. In mammals, the two light chains in an antibody are always identical, with only one type of light chain, kappa or lambda. The ratio of kappa to lambda is 70:30. However, with the occurrence of multiple myeloma or other B-cell malignancies this ratio is disturbed. Antibody to the kappa light chain is reportedly useful in the identification of leukemias, plasmacytomas, and certain non-Hodgkins lymphomas. Demonstration of clonality in lymphoid infiltrates indicates that the infiltrate is malignant.  |
| <b>Immunogen</b>               | Human B-Lymphoma Cells   |
| <b>Notes</b>                   | This product is provided under an intellectual property license from Life Technologies Corporation. The transfer of this product is conditioned on the buyer using the purchased product solely in research conducted by the buyer, excluding contract research or any fee for service research, and the buyer must not (1) use this product or its components for (a) diagnostic, therapeutic or prophylactic purposes; (b) testing, analysis or screening services, or information in return for compensation on a per-test basis; or (c) manufacturing or quality assurance or quality control, and/or (2) sell or transfer this product or its components for resale, whether or not resold for use in research. For information on purchasing a license to this product for purposes other than as described above, contact Life Technologies Corporation, 5781 Van Allen Way, Carlsbad, CA 92008 USA or <a href="mailto:outlicensing@thermofisher.com">outlicensing@thermofisher.com</a> . This conjugate is made on demand. Actual recovery may vary from the stated volume of this product. The volume will be greater than or equal to the unit size stated on the datasheet. |

**Product Application Details**

|                              |  |
|------------------------------|--|
| <b>Applications</b>          | Immunohistochemistry-Paraffin, Flow Cytometry, Immunohistochemistry, CyTOF-ready |
| <b>Recommended Dilutions</b> | Flow Cytometry, Immunohistochemistry, Immunohistochemistry-Paraffin, CyTOF-ready |
| <b>Application Notes</b>     | Optimal dilution of this antibody should be experimentally determined.           |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-34396AFP594**

---

|                    |  |
|--------------------|--|
| H00003514-P01-10ug | Recombinant Human Kappa Light Chain GST (N-Term) Protein |
| 243-B3-002         | TGF-beta 3 [Unconjugated]                                |
| NB7463             | Goat anti-Human Kappa Light Chain Secondary Antibody     |
| NB100-1028         | AIF-1/Iba1 Antibody                                      |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-34396AFP594](http://www.novusbio.com/reviews/submit/NBP2-34396AFP594)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

