## **Product Datasheet**

# E. Coli O157 Antibody (9.88) [DyLight 550] NBP2-54353R-0.1ml

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at: www.novusbio.com/NBP2-54353R

Updated 10/26/2023 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications
Submit a review at www.novusbio.com/reviews/destination/NBP2-54353R



## NBP2-54353R-0.1ml

F. Coli O157 Antibody (9.88) [Dyl ight 550].

E. Coli O157 Antibody (9.88) [DyLight 550]	
Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	9.88
Preservative	0.05% Sodium Azide
Isotype	IgM Kappa
Conjugate	DyLight 550
Purity	Protein L or PEG purified
Buffer	50mM Sodium Borate
<b>Product Description</b>	
Description	This conjugate is made on demand. Actual recovery may vary from the stated volume of this product. The volume will be greater than or equal to the unit size stated on the datasheet.
Host	Mouse
Species	E. coli
Reactivity Notes	E. coli O157
Specificity/Sensitivity	This monoclonal antibody shows specificity to E.coli O157 in a simple ELISA. Escherichia coli are Gram negative bacterium that are commonly found in the lower intestine of warm-blooded organisms (endotherms). Their serological types are determined in combination with somatic antigens (O group: O1-O173) and flagella antigens (H type: H1-H56). TheE. coli that cause intestinal infectious diseases including diarrhea, acute gastritis or colitis are referred to as pathogenicE. coli, which are classified into the following 4 groups according to differences in the mode of pathogenicity; enteropathogenicE. coli(EPEC), enteroinvasiveE. coli(EIEC), enterotoxigenicE. coli(ETEC) and enterohemorrhagicE. coli(EHEC). Although the identification of pathogenicE. coli requires verification of their pathogenicity, pathogenicE. coli often have specific serotypes; therefore, typing of the serogroup and serotype is necessary in screening pathogenicE. coli.
Immunogen	A crude sonicate of E. coli O157
Notes	DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.
<b>Product Application Details</b>	
Applications	ELISA, Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	ELISA, Immunocytochemistry/ Immunofluorescence

Product Application Details	
Applications	ELISA, Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	ELISA, Immunocytochemistry/ Immunofluorescence
Application Notes	Optimal dilution of this antibody should be experimentally determined.





### Novus Biologicals USA

10730 E. Briarwood Avenue Centennial, CO 80112 USA

Phone: 303.730.1950 Toll Free: 1.888.506.6887

Fax: 303.730.1966

nb-customerservice@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane Abingdon Science Park Abingdon, OX14 3NB, United Kingdom Phone: (44) (0) 1235 529449

Free Phone: 0800 37 34 15 Fax: (44) (0) 1235 533420 info.EMEA@bio-techne.com

#### **Bio-Techne Canada**

21 Canmotor Ave Toronto, ON M8Z 4E6 Canada

Phone: 905.827.6400 Toll Free: 855.668.8722 Fax: 905.827.6402

canada.inquires@bio-techne.com

### **General Contact Information**

www.novusbio.com Technical Support: nb-technical@biotechne.com

Orders: nb-customerservice@bio-techne.com

General: novus@novusbio.com

#### Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-54353R

Earn gift cards/discounts by submitting a publication using this product: www.novusbio.com/publications

