Product Datasheet

EGFR Antibody (h-R3 (Nimotuzumab)) [mFluor Violet 500 SE] NBP2-75186MFV500

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com

technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at: www.novusbio.com/NBP2-75186MFV500

Updated 9/20/2023 v.20.1

Earn rewards for product reviews and publications.



Submit a publication at www.novusbio.com/publications Submit a review at www.novusbio.com/reviews/destination/NBP2-75186MFV500

NBP2-75186MFV500

EGFR Antibody (h-R3 (Nimotuzumab)) [mFluor Violet 500 SE]

Product	Information
FIOUUCI	iniormation

Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	h-R3 (Nimotuzumab)
Preservative	0.05% Sodium Azide
Isotype	IgG1 Kappa
Conjugate	mFluor Violet 500 SE
Purity	Protein A purified
Buffer	50mM Sodium Borate
	•

Product Description

Human
1956
EGFR
Human
Binds to human EGFR domain III at a site distinct from antibody 7A7. (Talavera et al. 2011)
This antibody was generated by humanisation of the murine antibody R3, raised against EGFR purified from placenta membranes and A431 tumors.

Product Application Details

Applications	Flow Cytometry, Immunohistochemistry, Block/Neutralize, Immunofluorescence
Recommended Dilutions	Flow Cytometry, Immunohistochemistry, Immunofluorescence, Block/Neutralize
Application Notes	Optimal dilution of this antibody should be experimentally determined.

Images

EGFR Antibody (h-R3 (Nimotuzumab)) [mFluor Violet 500 SE] - Vial of mFluor Violet 500 conjugated antibody. mFluor Violet 500 is optimally excited at 410 nm by the Violet laser (405 nm) and has an emission maximum of 501 nm.

	mFluor™ \	/iolet 500	
	LASER (nm)	FILTER	
mFluor" Violet 500	Violet (405)	515/20	
	EXCITATION MAX (nm)	EMISSION MAX (nm)	Í
	410	501	
NOVUS		- 27-6331 	
CAUTION -Research Use Only			





Novus Biologicals USA

10730 E. Briarwood Avenue Centennial, CO 80112 USA Phone: 303.730.1950 Toll Free: 1.888.506.6887 Fax: 303.730.1966 nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave Toronto, ON M8Z 4E6 Canada Phone: 905.827.6400 Toll Free: 855.668.8722 Fax: 905.827.6402 canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane Abingdon Science Park Abingdon, OX14 3NB, United Kingdom Phone: (44) (0) 1235 529449 Free Phone: 0800 37 34 15 Fax: (44) (0) 1235 533420 info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com Technical Support: nb-technical@biotechne.com Orders: nb-customerservice@bio-techne.com General: novus@novusbio.com

Products Related to NBP2-75186MFV500

NBP3-06872MFV500Human IgG1 Kappa Isotype Control [mFluor Violet 500 SE]NBP1-84815PEPEGFR Recombinant Protein Antigen236-EG-200EGF [Unconjugated]DYC1095B-2EGFR [Biotin]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-75186MFV500

Earn gift cards/discounts by submitting a publication using this product: www.novusbio.com/publications

www.novusbio.com

