

Product Datasheet

Kidney Matched Tumor and Normal Tissue Slides (Tumor) NBP2-77671

Unit Size: 10 Slides

Store at -70C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-77671

Updated 3/5/2020 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-77671




NBP2-77671**Kidney Matched Tumor and Normal Tissue Slides (Tumor)**

Product Information	
Unit Size	10 Slides
Concentration	Concentration is not relevant for this product. Please see the protocols for proper use of this product.
Storage	Store at -70C. Avoid freeze-thaw cycles.
Product Description	
Description	Frozen tissue slides are ideal for rapid cellular localization of RNA and protein expression, especially for antibodies that don't recognize paraffin embedded tissue. Tissues were snap frozen in liquid nitrogen immediately after excision, and embedded in OCT. Frozen tissue sections in each pair are prepared from the same donor. Tissue sections are 5-10 m in thickness, mounted on positively charged glass slides, and fixed by cold acetone. At least one of the tissue slides from each lot was stained with H & E for quality assurance. Product is supplied as 5 diseased slides and 5 adjacent normal slides.
Species	Human
Notes	Donor information available upon request
Lysate Type	Matched Tumor and Normal Tissue
Lysate Tissue	Kidney
Lysate Tissue Condition	Tumor
Product Application Details	
Applications	Hematoxylin and Eosin Stain, Immunohistochemistry, Immunohistochemistry-Frozen, In-situ Hybridization
Recommended Dilutions	Immunohistochemistry, Immunohistochemistry-Frozen, In-situ Hybridization, Hematoxylin and Eosin Stain
Application Notes	Each individual tissue spot should be >90% core retention and total core lost from 1 slide should be no more than 5%




Images

Hematoxylin & Eosin Stain: Kidney Matched Tumor & Normal Tissue Slides (Tumor) [NBP2-77671] - Tissue: Human Kidney Tumor, Pathology: Primary Tumor



Hematoxylin & Eosin Stain: Kidney Matched Tumor & Normal Tissue Slides (Tumor) [NBP2-77671] - Tissue: Human Kidney, Pathology: Adjacent normal tissue





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-77671

NBP2-30220	Human Kidney Tissue MicroArray (Cancer)
------------	---

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Slides are guaranteed for 3 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-77671

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

