# **Product Datasheet**

# Abdominal Tissue Slides (Neurofibromatosis)- Paraffin NBP2-77854

Unit Size: 5 Slides

Store at 4C.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at: www.novusbio.com/NBP2-77854

Updated 1/9/2024 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications Submit a review at www.novusbio.com/reviews/destination/NBP2-77854



#### NBP2-77854

Abdominal Tissue Slides (Neurofibromatosis)- Paraffin


Abdominal Floods Charles (Notice Indiana)		
Product Information		
Unit Size	5 Slides	
Concentration	Concentration is not relevant for this product. Please see the protocols for proper use of this product.	
Storage	Store at 4C.	
<b>Product Description</b>		

Storage	Store at 4C.
<b>Product Description</b>	
Description	FFPE tissue slides are ideal for rapid cellular localization of RNA and protein expression.  Tissues were fixed in formalin immediately after excision, and embedded in paraffin. Tissue sections are 5um in thickness, and mounted on positively charged glass slides. At least one of the tissue slides from each lot was stained with H & E for quality assurance.
Species	Human
Notes	Donor information available upon request
Lysate Type	Tissue
Lysate Tissue	Abdominal
Lysate Tissue Condition	Neurofibromatosis

_	
<b>Product Application Details</b>	
Applications	Hematoxylin and Eosin Stain, Immunohistochemistry, Immunohistochemistry- Paraffin, In-situ Hybridization
Recommended Dilutions	Immunohistochemistry, Immunohistochemistry-Paraffin, In-situ Hybridization, Hematoxylin and Eosin Stain
Application Notes	Each individual tissue spot should be >90% core retention and total core lost from 1 slide should be no more than 5%. Please bake slides at 60C for 30 minutes before use.

# **Images**

Hematoxylin & Eosin Stain: Abdominal Tissue Slides (Neurofibromatosis) [NBP2-77854] - Tissue: Human Abdomen, Pathology: Neurofibromatosis







# Novus Biologicals USA

10730 E. Briarwood Avenue Centennial, CO 80112 USA

Phone: 303.730.1950 Toll Free: 1.888.506.6887

Fax: 303.730.1966

nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave Toronto, ON M8Z 4E6 Canada

Phone: 905.827.6400 Toll Free: 855.668.8722 Fax: 905.827.6402

canada.inquires@bio-techne.com

#### **Bio-Techne Ltd**

19 Barton Lane Abingdon Science Park Abingdon, OX14 3NB, United Kingdom Phone: (44) (0) 1235 529449

Free Phone: 0800 37 34 15 Fax: (44) (0) 1235 533420 info.EMEA@bio-techne.com

#### **General Contact Information**

www.novusbio.com Technical Support: nb-technical@biotechne.com

Orders: nb-customerservice@bio-techne.com

General: novus@novusbio.com

#### Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Slides are guaranteed for 3 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-77854

Earn gift cards/discounts by submitting a publication using this product: www.novusbio.com/publications

