Product Datasheet

Pelvic Tissue Slides (Malignant Mesothelioma)- Paraffin NBP2-77958

Unit Size: 5 Slides

Store at 4C.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at: www.novusbio.com/NBP2-77958

Updated 1/9/2024 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications Submit a review at www.novusbio.com/reviews/destination/NBP2-77958



NBP2-77958

Pelvic Tissue Slides (Malignant Mesothelioma)- Paraffin


Charles (manighter and charles)		
Product Information		
Unit Size	5 Slides	
Concentration	Concentration is not relevant for this product. Please see the protocols for proper use of this product.	
Storage	Store at 4C.	
Product Description		

Storage	Store at 4C.	
Product Description		
Description	FFPE tissue slides are ideal for rapid cellular localization of RNA and protein expression. Tissues were fixed in formalin immediately after excision, and embedded in paraffin. Tissue sections are 5um in thickness, and mounted on positively	
	charged glass slides. At least one of the tissue slides from each lot was stained with H & E for quality assurance.	
Species	Human	
Notes	Donor information available upon request	
Lysate Type	Tissue	
Lysate Tissue	Pelvic	
Lysate Tissue Condition	Malignant Mesothelioma	

Product Application Details	
Applications	Hematoxylin and Eosin Stain, Immunohistochemistry, Immunohistochemistry-Paraffin, In-situ Hybridization
Recommended Dilutions	Immunohistochemistry, Immunohistochemistry-Paraffin, In-situ Hybridization, Hematoxylin and Eosin Stain
Application Notes	Each individual tissue spot should be >90% core retention and total core lost from 1 slide should be no more than 5%. Please bake slides at 60C for 30 minutes before use.

Images

Hematoxylin & Eosin Stain: Pelvic Tissue Slides (Malignant Mesothelioma) [NBP2-77958] - Tissue: Human Pelvis, Pathology: Malignant Mesothelioma





Novus Biologicals USA

10730 E. Briarwood Avenue Centennial, CO 80112 USA

Phone: 303.730.1950 Toll Free: 1.888.506.6887

Fax: 303.730.1966

nb-customerservice@bio-techne.com

Bio-Techne Ltd

19 Barton Lane Abingdon Science Park Abingdon, OX14 3NB, United Kingdom Phone: (44) (0) 1235 529449

Free Phone: 0800 37 34 15 Fax: (44) (0) 1235 533420 info.EMEA@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave Toronto, ON M8Z 4E6 Canada

Phone: 905.827.6400 Toll Free: 855.668.8722 Fax: 905.827.6402

canada.inquires@bio-techne.com

General Contact Information

www.novusbio.com Technical Support: nb-technical@biotechne.com

Orders: nb-customerservice@bio-techne.com

General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Slides are guaranteed for 3 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-77958

Earn gift cards/discounts by submitting a publication using this product: www.novusbio.com/publications

