# **Product Datasheet**

# Human FDA Standard Tissue MicroArray (Normal) NBP2-78057

Unit Size: 1 Slides

Store at 4C.

www.novusbio.com



technical@novusbio.com

**Publications: 1** 

Protocols, Publications, Related Products, Reviews, Research Tools and Images at: www.novusbio.com/NBP2-78057

Updated 3/7/2024 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications Submit a review at www.novusbio.com/reviews/destination/NBP2-78057



# NBP2-78057

Human FDA Standard Tissue MicroArray (Normal)

Product Information	
Unit Size	1 Slides
Concentration	Concentration is not relevant for this product. Please see the protocols for proper use of this product.
Storage	Store at 4C.
Product Description	
Description	Please see online datasheet for well details: <a href="https://www.novusbio.com/NBP2-78057">www.novusbio.com/NBP2-78057</a>
Species	Human
Lysate Type	Tissue
Lysate Tissue	Multi-tissue
Lysate Tissue Condition	Normal
Product Application Details	
Applications	Immunohistochemistry, Immunohistochemistry-Paraffin, In-situ Hybridization, RNAscope compatible
Recommended Dilutions	Immunohistochemistry, Immunohistochemistry-Paraffin, In-situ Hybridization, RNAscope compatible
Application Notes	Each individual tissue spot should be >90% core retention and total core lost from 1 slide should be no more than 5%. Please bake slides at 60C for 30 minutes before use.

Complete product information can be found online at www.novusbio.com/NBP2-78057. Please contact technical service with any questions regarding the use of this product.

#### **Publications**

Hasegawa M, Hanamatsu Y, Saigo C et al. Tumor suppressor effect of an antibody on xenotransplanted sarcomatoid mesothelioma cells Thoracic cancer 2022-08-02 [PMID: 35916358] (IF/IHC, Human)





# **Novus Biologicals USA**

10730 E. Briarwood Avenue Centennial, CO 80112 USA

Phone: 303.730.1950 Toll Free: 1.888.506.6887

Fax: 303.730.1966

nb-customerservice@bio-techne.com

# Bio-Techne Canada 21 Canmotor Ave

21 Canmotor Ave Toronto, ON M8Z 4E6 Canada

Phone: 905.827.6400 Toll Free: 855.668.8722 Fax: 905.827.6402

canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane Abingdon Science Park Abingdon, OX14 3NB, United Kingdom Phone: (44) (0) 1235 529449

Free Phone: 0800 37 34 15 Fax: (44) (0) 1235 533420 info.EMEA@bio-techne.com

# **General Contact Information**

www.novusbio.com Technical Support: nb-technical@biotechne.com

Orders: nb-customerservice@bio-techne.com

General: novus@novusbio.com

#### Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Tissue Micro Arrays are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-78057

Earn gift cards/discounts by submitting a publication using this product: www.novusbio.com/publications

