Product Datasheet

Human Multiple Tumor Tissue MicroArray (Tumor) NBP2-78109

Unit Size: 5 Slides

Store at 4C.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at: www.novusbio.com/NBP2-78109

Updated 5/8/2019 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications Submit a review at www.novusbio.com/reviews/destination/NBP2-78109



NBP2-78109

Human Multiple Tumor Tissue MicroArray (Tumor)

Product Information	
5 Slides	
Concentration is not relevant for this product. Please see the protocols for proper use of this product.	
Store at 4C.	
Product Description	
Please see online datasheet for well details: www.novusbio.com/NBP2-78109	
Human	
Tissue	
Multi-tissue	
Tumor	
Immunohistochemistry, Immunohistochemistry-Paraffin, In-situ Hybridization	
Immunohistochemistry, Immunohistochemistry-Paraffin, In-situ Hybridization	
Each individual tissue spot should be >90% core retention and total core lost from 1 slide should be no more than 5%. Please bake slides at 60C for 30 minutes before use.	

Complete product information can be found online at www.novusbio.com/NBP2-78109. Please contact technical service with any questions regarding the use of this product.





Novus Biologicals USA

10730 E. Briarwood Avenue Centennial, CO 80112 USA

Phone: 303.730.1950 Toll Free: 1.888.506.6887

Fax: 303.730.1966

nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave Toronto, ON M8Z 4E6 Canada

Phone: 905.827.6400 Toll Free: 855.668.8722 Fax: 905.827.6402

canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane Abingdon Science Park Abingdon, OX14 3NB, United Kingdom Phone: (44) (0) 1235 529449

Free Phone: 0800 37 34 15 Fax: (44) (0) 1235 533420 info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com Technical Support: nb-technical@biotechne.com

Orders: nb-customerservice@bio-techne.com

General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Tissue Micro Arrays are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-78109

Earn gift cards/discounts by submitting a publication using this product: www.novusbio.com/publications

