Product Datasheet

Fluorescein Antibody (4-4-20 (enhanced)) [mFluor Violet 500 SE] NBP2-81026MFV500

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at: www.novusbio.com/NBP2-81026MFV500

Updated 9/20/2023 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications
Submit a review at www.novusbio.com/reviews/destination/NBP2-81026MFV500



NBP2-81026MFV500

Fluorescein Antibody (4-4-20 (enhanced)) [mFluor Violet 500 SE]

| · ···································· | |
|--|---|
| Product Information | |
| Unit Size | 0.1 ml |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. |
| Clonality | Monoclonal |
| Clone | 4-4-20 (enhanced) |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG1 Kappa |
| Conjugate | mFluor Violet 500 SE |
| Purity | Protein A purified |
| Buffer | 50mM Sodium Borate |
| Product Description | |
| Host | Rat |
| Species | Non-species specific |
| Specificity/Sensitivity | Fluorescein labelled proteins. |
| Immunogen | Fluorescein I-bovine serum albumin (FII-BSA). |
| Product Application Details | |
| Applications | Western Blot, Flow Cytometry, Control |
| Recommended Dilutions | Western Blot, Flow Cytometry, Control |

Optimal dilution of this antibody should be experimentally determined.

Images

Application Notes

Fluorescein Antibody (4-4-20 (enhanced)) [mFluor Violet 500 SE] - Vial of mFluor Violet 500 conjugated antibody. mFluor Violet 500 is optimally excited at 410 nm by the Violet laser (405 nm) and has an emission maximum of 501 nm.





Novus Biologicals USA

10730 E. Briarwood Avenue Centennial, CO 80112 USA

Phone: 303.730.1950 Toll Free: 1.888.506.6887

Fax: 303.730.1966

nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave Toronto, ON M8Z 4E6

Canada

Phone: 905.827.6400 Toll Free: 855.668.8722 Fax: 905.827.6402

canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane Abingdon Science Park Abingdon, OX14 3NB, United Kingdom Phone: (44) (0) 1235 529449

Free Phone: 0800 37 34 15 Fax: (44) (0) 1235 533420 info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com

Technical Support: nb-technical@bio-

techne.com

Orders: nb-customerservice@bio-techne.com

General: novus@novusbio.com

Products Related to NBP2-81026MFV500

NBP2-33376H Blue Marker Antibody (6F4-F6) [HRP]

NBP1-51817MFV500 Rat IgG1 Kappa Light Chain Isotype Control (KLH/G1-2-2) [mFluor

Violet 500 SE]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-81026MFV500

Earn gift cards/discounts by submitting a publication using this product: www.novusbio.com/publications

