# **Product Datasheet**

## EphA2 Antibody (1C1) [mFluor Violet 610 SE] NBP2-81069MFV610

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com

technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at: www.novusbio.com/NBP2-81069MFV610

Updated 9/20/2023 v.20.1







#### NBP2-81069MFV610

EphA2 Antibody (1C1) [mFluor Violet 610 SE]

Product Information			
Unit Size	0.1 ml		
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.		
Storage	Store at 4C in the dark.		
Clonality	Monoclonal		
Clone	1C1		
Preservative	0.05% Sodium Azide		
Isotype	IgG1 Kappa		
Conjugate	mFluor Violet 610 SE		
Purity	Protein A purified		
Buffer	50mM Sodium Borate		
Product Description			
Host	Human		
Gene ID	1969		
Gene Symbol	EPHA2		
Species	Human		
Specificity/Sensitivity	1C1 binds specifically to the amino-termianal LBD (ligand-binding domain) region of the EphA2 ectodomain (IgG1 KD ~0.8 nM; Fab KD ~140 nM) and not to other related human Eph receptors such as EphA1, EphA2, and EphA4, or any EphB receptor family members.		
Immunogen	1C1 was isolated after screening of a non-immunized human Fab phage library. Selection was based on the ability to bind soluble human recombinant EphA2 (residues 183-510). The Fab was subsequently converted to full length human IgG1.		
Product Application Details			
Applications	Western Blot, ELISA, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Block/Neutralize		
Recommended Dilutions	Western Blot, Flow Cytometry, ELISA, Immunocytochemistry/ Immunofluorescence, Block/Neutralize		

### Images

**Application Notes** 

EphA2 Antibody (1C1) [mFluor Violet 610 SE] [NBP2-81069MFV610] - Vial of mFluor Violet 610 conjugated antibody. mFluor Violet 610 is optimally excited at 421 nm by the Violet laser (405 nm) and has an emission maximum of 613 nm.

	mFluor™∖	/iolet 610	
and the second s	LASER (nm)	FILTER	1
<sup>mFluor<sup>™</sup> Violet 6<sup>10</sup></sup>	Violet (405)	605/30	
	EXCITATION MAX (nm)	EMISSION MAX (nm)	1
<b>G</b> NOVUS	421	613	

www.novusbio.com



Optimal dilution of this antibody should be experimentally determined.



#### **Novus Biologicals USA**

10730 E. Briarwood Avenue Centennial, CO 80112 USA Phone: 303.730.1950 Toll Free: 1.888.506.6887 Fax: 303.730.1966 nb-customerservice@bio-techne.com

#### **Bio-Techne Canada**

21 Canmotor Ave Toronto, ON M8Z 4E6 Canada Phone: 905.827.6400 Toll Free: 855.668.8722 Fax: 905.827.6402 canada.inquires@bio-techne.com

#### **Bio-Techne Ltd**

19 Barton Lane Abingdon Science Park Abingdon, OX14 3NB, United Kingdom Phone: (44) (0) 1235 529449 Free Phone: 0800 37 34 15 Fax: (44) (0) 1235 533420 info.EMEA@bio-techne.com

#### **General Contact Information**

www.novusbio.com Technical Support: nb-technical@biotechne.com Orders: nb-customerservice@bio-techne.com General: novus@novusbio.com

#### Products Related to NBP2-81069MFV610

639-A2-200	EphA2 [Unconjugated]
DVE00	VEGF [HRP]
NBP3-06872MFV610	Human IgG1 Kappa Isotype Control [mFluor Violet 610 SE]
NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]

#### Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-81069MFV610

Earn gift cards/discounts by submitting a publication using this product: www.novusbio.com/publications

www.novusbio.com



www.novusbio.com



technical@novusbio.com