Product Datasheet

TIGIT Antibody (4D4) [mFluor Violet 450 SE] NBP2-81116MFV450

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at: www.novusbio.com/NBP2-81116MFV450

Updated 5/1/2025 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications
Submit a review at www.novusbio.com/reviews/destination/NBP2-81116MFV450



NBP2-81116MFV450

T Antibody (4D4) [mFluor Violet 450 SF]

| TIGIT Antibody (4D4) [mFluor Violet 450 SE] | |
|---|---|
| Product Information | |
| Unit Size | 0.1 ml |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. |
| Clonality | Monoclonal |
| Clone | 4D4 |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG Kappa |
| Conjugate | mFluor Violet 450 SE |
| Purity | Protein A purified |
| Buffer | 50mM Sodium Borate |
| Product Description | |
| Host | Armenian Hamster |
| Gene ID | 201633 |
| Gene Symbol | TIGIT |
| Species | Mouse |
| Specificity/Sensitivity | 4D4 specifically recognises TIGIT. |
| Immunogen | Armenian hamsters were immunized with recombinant mouse TIGIT tetramers by a combination of s.c. and food pad immunization and booster injections. Draining lymph nodes were fused with Sp2/0-Ag14 cells, selected in HAT (hypoxanthine/aminopterin/thymidine) medium and supernatants were screened for specificity by ELISA and flow cytometry using TIGIT-transfectants. |
| Product Application Details | |
| Applications | Flow Cytometry, Block/Neutralize |
| Recommended Dilutions | Flow Cytometry, Block/Neutralize |
| Annilostian Natas | Onting all dilutions of this partitional to be a labelle and a supering and the determined |

Application Notes Optimal dilution of this antibody should be experimentally determined.

Images

TIGIT Antibody (4D4) [mFluor Violet 450 SE] [NBP2-81116MFV450] -Vial of mFluor Violet 450 conjugated antibody. mFluor Violet 450 is optimally excited at 406 nm by the Violet laser (405 nm) and has an emission maximum of 445 nm.







Novus Biologicals USA

10730 E. Briarwood Avenue Centennial, CO 80112

USA

Phone: 303.730.1950 Toll Free: 1.888.506.6887

Fax: 303.730.1966

nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave Toronto, ON M8Z 4E6

Canada

Phone: 905.827.6400 Toll Free: 855.668.8722 Fax: 905.827.6402

canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane Abingdon Science Park Abingdon, OX14 3NB, United Kingdom

Phone: (44) (0) 1235 529449 Free Phone: 0800 37 34 15 Fax: (44) (0) 1235 533420 info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com

Technical Support: nb-technical@bio-

techne.com

Orders: nb-customerservice@bio-techne.com

General: novus@novusbio.com

Products Related to NBP2-81116MFV450

DY417-05 IL-10 [Biotin]

7898-TGB-100 TIGIT [Unconjugated]

AF1086 PD-1 Antibody [Unconjugated]
NBP1-19371 CD4 Antibody - BSA Free

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-81116MFV450

Earn gift cards/discounts by submitting a publication using this product: www.novusbio.com/publications

