

Product Datasheet

desmosome/hemidesmosome Antibody (F12) - Azide and BSA Free NBP2-81279-0.05mg

Unit Size: 0.05 mg

Store at 4C for up to 3 months. For longer storage, aliquot and store at -20C.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-81279

Updated 4/19/2023 v.20.1

**Earn rewards for product
reviews and publications.**

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-81279



NBP2-81279-0.05mg

desmosome/hemidesmosome Antibody (F12) - Azide and BSA Free

Product Information	
Unit Size	0.05 mg
Concentration	1 mg/ml
Storage	Store at 4C for up to 3 months. For longer storage, aliquot and store at -20C.
Clonality	Monoclonal
Clone	F12
Preservative	0.02% Proclin 300
Isotype	IgM Kappa
Purity	Affinity purified
Buffer	PBS

Product Description	
Host	Human
Species	Human
Specificity/Sensitivity	"F12 recognizes an unknown polypeptide of desmosomal and hemidesmosomal plaques with a MW of 185 kDa (determined from immunoblotting using bovine tongue epithelium extract as the substrate), and by IEM (immunoelectron microscopy) it was determined that it is the intracellular part of the desmosome that is bound (Gilbert, 1997). F12 shares some immunochemical properties with autoantibodies present in paraneoplastic pemphigus sera. Pemphigus diseases are those resulting in blistering of the skin due to the production of autoantibodies against keratinocyte surface proteins - autoantibodies of patients with pemphigus vulgaris recognize desmoglein 1 (Dsg1), a major glycoprotein of the desmosome complex."
Immunogen	F12 was established by fusing peripheral blood lymphocytes derived from a human vulgaris patient with a heteromyeloma cell line.

Product Application Details	
Applications	Electron Microscopy, Flow Cytometry, Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Electron Microscopy
Application Notes	F12 labels the keratinocyte membrane and the basement membrane zone of rat tongue sections by indirect immunofluorescence (IIF). The antibody also stains both the cell membrane and the basement membrane zone of stratified squamous epithelia, and other epithelial tissues such as urinary bladder, small bowel, thymus, and liver. Non-epithelial tissues, such as myocardium are additionally stained. Sera from patients with pemphigus vulgaris can be incubated with human skin sections to determine their capacity to block the binding of F12 in an IIF blocking assay. IEM (immunoelectron microscopy) can be used to determine the ultralocalization of the antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-81279-0.05mg

G-102-C	Goat anti-Human IgG Secondary Antibody [Unconjugated]
NB7446	Goat anti-Human IgG Fc Secondary Antibody

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-81279

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

