Product Datasheet

MST1/STK4 Antibody (01) [mFluor Violet 450 SE] NBP3-06376MFV450

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com

technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at: www.novusbio.com/NBP3-06376MFV450

Updated 7/11/2023 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications Submit a review at www.novusbio.com/reviews/destination/NBP3-06376MFV450



NBP3-06376MFV450

MST1/STK4 Antibody (01) [mFluor Violet 450 SE]

,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,			
Product Information			
Unit Size	0.1 ml		
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.		
Storage	Store at 4C in the dark.		
Clonality	Monoclonal		
Clone	01		
Preservative	0.05% Sodium Azide		
Isotype	IgG1		
Conjugate	mFluor Violet 450 SE		
Purity	Protein A purified		
Buffer	50mM Sodium Borate		
Product Description			
Host	Mouse		
Species	Human		
Immunogen	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human MST1/STK4 (Uniprot#: Q13043-1; Glu2-Phe487).		
Notes	mFluor(TM) is a trademark of AAT Bioquest, Inc. This conjugate is made on demand. Actual recovery may vary from the stated volume of this product. The volume will be greater than or equal to the unit size stated on the datasheet.		
Declarat Annihestican Detaile			

Product Application Details

Applications	Flow Cytometry, Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	Flow Cytometry, Immunocytochemistry/ Immunofluorescence
Application Notes	Optimal dilution of this antibody should be experimentally determined.

Images

MST1/STK4 Antibody (01) [mFluor Violet 450 SE] [NBP3-06376MFV450] - Vial of mFluor Violet 450 conjugated antibody. mFluor Violet 450 is optimally excited at 406 nm by the Violet laser (405 nm) and has an emission maximum of 445 nm.

The second second				
Contraction of the second second second				
	mFluor™∖	/iolet 450		
Same and the second sec			1	
	LASER (nm)	FILTER	-	
^I ^{™Fluor™} Violet 450	Violet (405)	450/45		
			1	
	EXCITATION MAX (nm)	EMISSION MAX (nm)		
	406	445		
biotechne brand				
CAUTION -Research Use Dri				





Novus Biologicals USA

10730 E. Briarwood Avenue Centennial, CO 80112 USA Phone: 303.730.1950 Toll Free: 1.888.506.6887 Fax: 303.730.1966 nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave Toronto, ON M8Z 4E6 Canada Phone: 905.827.6400 Toll Free: 855.668.8722 Fax: 905.827.6402 canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane Abingdon Science Park Abingdon, OX14 3NB, United Kingdom Phone: (44) (0) 1235 529449 Free Phone: 0800 37 34 15 Fax: (44) (0) 1235 533420 info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com Technical Support: nb-technical@biotechne.com Orders: nb-customerservice@bio-techne.com General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-06376MFV450

Earn gift cards/discounts by submitting a publication using this product: www.novusbio.com/publications

www.novusbio.com

