

# Product Datasheet

## CD9 Immunobeads for Exosome Isolation NBP3-11768

Unit Size: 10 Reactions

Store at 4C.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-11768](http://www.novusbio.com/NBP3-11768)

Updated 7/24/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-11768](http://www.novusbio.com/reviews/destination/NBP3-11768)



**NBP3-11768****CD9 Immunobeads for Exosome Isolation**

<b>Product Information</b>	
<b>Unit Size</b>	10 Reactions
<b>Concentration</b>	Please see the protocols for proper use of this product. If no protocol is available, contact technical services for assistance.
<b>Storage</b>	Store at 4C.
<b>Product Description</b>	
<b>Description</b>	CD9 Immunocapture Beads are 0.4 um latex beads covalently coupled with antibodies against extracellular vesicle (EV) surface antigen CD9. Samples containing EVs are incubated with Immunocapture Beads for rapid and specific isolation of EVs expressing CD9. Using the supplied Elution Buffer, EVs can be eluted from Immunocapture beads and used in downstream analyses such as protein profiling and nucleic acid extraction. Immunobeads can be re-used after incubation with Bead Regeneration Buffer.
<b>Species</b>	Human
<b>Kit Components</b>	CD9 Immunocapture Beads Exosome Elution Buffer Bead Regeneration Buffer
<b>Suitable Sample Type</b>	Biofluids (serum, plasma, and urine)
<b>Sample Volume</b>	0.5-1 mL





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Support products are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-11768](http://www.novusbio.com/reviews/submit/NBP3-11768)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

