

# Product Datasheet

## Donor Bovine Serum S11350-NOV

Unit Size: 500 ml

Store at -20C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/S11350-NOV](http://www.novusbio.com/S11350-NOV)

Updated 6/1/2021 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/S11350-NOV](http://www.novusbio.com/reviews/destination/S11350-NOV)



**S11350-NOV**

## Donor Bovine Serum

**Product Information**

<b>Unit Size</b>	500 ml
<b>Concentration</b>	Please see the protocols for proper use of this product. If no protocol is available, contact technical services for assistance.
<b>Storage</b>	Store at -20C. Avoid freeze-thaw cycles.

**Product Description**

<b>Description</b>	<p>Donor Bovine Serum Summary</p> <p>Donor Bovine Serum (DBS) is frequently used in applications where serum sero-negative to IBR and BVD is required.</p> <p>Why Use Donor Bovine Serum</p> <p>Donor Bovine Serum is collected from healthy, live animals greater 12 months of age or older. The donor cattle are carefully selected and maintained in closed herds ensuring complete control of their care, diet and health. Before acceptance into the segregated herd, all cattle must test sero-negative to Infectious Bovine Rhinotracheitis (IBR) and Bovine Viral Diarrhea (BVD). The animals are used exclusively to produce donor bovine products such as bovine serum and are not medicated or fed sub-therapeutic levels of antibiotics.</p> <p>Specifications</p> <p>Purification: 0.2 um filtered</p> <p>Sterility Testing: Pass</p> <p>Country of Origin: UNITED STATES</p> <p>Limitations</p> <p>For research use only. Not for diagnostic use.</p>
<b>Species</b>	Bovine

**Product Application Details**

<b>Applications</b>	Cell Culture
<b>Recommended Dilutions</b>	Cell Culture
<b>Application Notes</b>	Donor Bovine Serum (DBS) is frequently used in applications where serum sero-negative to IBR and BVD is required.



## Images

Serum: Donor Bovine Serum [S11350-NOV] - Serum: Donor Bovine Serum [s11350-NOV] - Donor Bovine Serum





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-  
techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Support products are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/S11350-NOV](http://www.novusbio.com/reviews/submit/S11350-NOV)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

