

# Product Datasheet

## Tripeptidyl-Peptidase I/TPP1 Antibody H00001200-B01

Unit Size: 0.05 ml

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/H00001200-B01](http://www.novusbio.com/H00001200-B01)

Updated 2/5/2017 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/H00001200-B01](http://www.novusbio.com/reviews/destination/H00001200-B01)



**H00001200-B01****Tripeptidyl-Peptidase I/TPP1 Antibody**

<b>Product Information</b>	
<b>Unit Size</b>	0.05 ml
<b>Concentration</b>	This product is unpurified. The exact concentration of antibody is not quantifiable.
<b>Storage</b>	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
<b>Clonality</b>	Polyclonal
<b>Preservative</b>	No Preservative
<b>Purity</b>	Unpurified
<b>Buffer</b>	This is unpurified antiserum so it is not possible to supply the specific antibody concentration.

<b>Product Description</b>	
<b>Host</b>	Mouse
<b>Gene ID</b>	1200
<b>Gene Symbol</b>	TPP1
<b>Species</b>	Human
<b>Reactivity Notes</b>	Human. Other species not tested.
<b>Specificity/Sensitivity</b>	TPP1 - tripeptidyl peptidase I,
<b>Immunogen</b>	TPP1 (NP_000382, 1 a.a. ~ 563 a.a) full-length human protein.MGLQACLLGLFALILSGKCSYSPEPDQRRTLPPGWVSLGRADPEEELS LTFALRQQNVERLSELVQAVSDPSSPQYGKYLLENVADLVRPSPLTLHTVQK WLLAAGAQQKCHSVITQDFLTCWLSIRQAELLPLGAEFHYYVGGPTETHVVRSP HPYQLPQALAPHVDFVGGGLHRFPPTSSLRQRPEPQVTGTVGLHLGVTPSVIRK RYNLTSQDVGSGTSNNSQACAQFLEQYFHDSDLAQFMRLFGGNFAHQA
<b>Notes</b>	This product is produced by and distributed for Abnova, a company based in Taiwan.

<b>Product Application Details</b>	
<b>Applications</b>	Western Blot, ELISA, Immunohistochemistry, Immunohistochemistry-Paraffin
<b>Recommended Dilutions</b>	Western Blot 1:500, ELISA, Immunohistochemistry, Immunohistochemistry-Paraffin
<b>Application Notes</b>	Antibody reactive against transfected lysate and tissue lysate for Western Blot. Has also been used for immunohistochemistry (paraffin) and ELISA.



### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/H00001200-B01](http://www.novusbio.com/reviews/submit/H00001200-B01)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

