

Product Datasheet

Recombinant Human Kininogen GST (N-Term) Protein H00003827-Q01-25ug

Unit Size: 25 ug

Store at -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00003827-Q01

Updated 10/18/2022 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00003827-Q01




H00003827-Q01-25ug**Recombinant Human Kininogen GST (N-Term) Protein**

| Product Information | |
|------------------------------------|--|
| Unit Size | 25 ug |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at -80C. Avoid freeze-thaw cycles. |
| Preservative | No Preservative |
| Purity | >80% by SDS-PAGE and Coomassie blue staining |
| Buffer | 50 mM Tris-HCl, 10 mM reduced Glutathione, pH 8.0 in the elution buffer. |
| Target Molecular Weight | 37.4 kDa |
| Product Description | |
| Description | <p>A recombinant protein with a N-terminal GST tag corresponding to the amino acid sequence 322-427 of Human Kininogen</p> <p>Source: <i>Wheat Germ (in vitro)</i></p> <p>Amino Acid Sequence: VARETTCSKESNEELTESCETKKLGGQSLDCNAEVYVVPWEKKIYPTVNCQPLG MISLMKRPPGFPFRSSRIGEIKEETTSHLRSCEYKGRPPKAGAEPASEREVS</p> |
| Gene ID | 3827 |
| Gene Symbol | KNG1 |
| Species | Human |
| Preparation Method | in vitro wheat germ expression system |
| Details of Functionality | This protein was produced in an in vitro wheat germ expression system that should preserve correct conformational folding that is necessary for biological function. While it is possible that this protein could display some level of activity, the functionality of this protein has not been explicitly measured or validated. |
| Notes | This product is produced by and distributed for Abnova, a company based in Taiwan. |
| Product Application Details | |
| Applications | Western Blot, ELISA, Protein Array, SDS-Page, Immunoaffinity Purification |
| Recommended Dilutions | Western Blot, ELISA, SDS-Page, Protein Array, Immunoaffinity Purification |

Images

SDS-Page: Recombinant Human Kininogen GST (N-Term) Protein [H00003827-Q01] - 12.5% SDS-PAGE Stained with Coomassie Blue.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. This product is guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00003827-Q01

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

