

Product Datasheet

Apolipoprotein F Antibody

H00000319-B01P

Unit Size: 0.05 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00000319-B01P

Updated 12/26/2016 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00000319-B01P



H00000319-B01P

Apolipoprotein F Antibody


| Product Information | |
|---------------------|--|
| Unit Size | 0.05 mg |
| Concentration | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| Storage | Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles. |
| Clonality | Polyclonal |
| Preservative | No Preservative |
| Purity | Protein G purified |
| Buffer | PBS (pH 7.4) |

| Product Description | |
|---------------------|--|
| Host | Mouse |
| Gene ID | 319 |
| Gene Symbol | APOF |
| Species | Human |
| Immunogen | APOF (AAH26257, 37 a.a. - 327 a.a.) full-length human protein. TSYGKQTNVLMHFPLSLESQTPSSDPLSCQFLHPKSLPGFSHMAPLPKFLVSL ALRNDALEEAGCQADVWALQLQLYRQGGVNATQVLIQHRLRGLQKGRSTERNV SVEALASALQLLAREQQSTGRVGRSLPTEDCENEKEQAVHNVVQLLPGVGTFY NLGTALYYATQNCLGKARERGRDGAIDLGYDLLMTMAGMSGGPMGLAISAALK PALRSGVQQLIQYYQDQKDANISQPETTKEGLRAISDVSDLEETTTLASFISEVV SSAPYWGWAIIKSYDLDPGAGSLEI |
| Notes | Quality control test: Antibody reactive against mammalian transfected lysate. This product is produced by and distributed for Abnova, a company based in Taiwan. |

| Product Application Details | |
|-----------------------------|---|
| Applications | Western Blot |
| Recommended Dilutions | Western Blot |
| Application Notes | It has been used for WB and Functional. |

Images

Western Blot: Apolipoprotein F Antibody [H00000319-B01P] - Analysis of APOF expression in transfected 293T cell line by APOF polyclonal antibody. Lane 1: APOF transfected lysate(32.01 KDa). Lane 2: Non-transfected lysate.





Novus Biologicals USA

8100 Southpark Way, A-8
Littleton, CO 80120
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Novus Biologicals Canada

461 North Service Road West, Unit B37
Oakville, ON L6M 2V5
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada@novusbio.com

Novus Biologicals Europe

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Products Related to H0000319-B01P

| | |
|----------------------|---|
| HAF007 | Goat anti-Mouse IgG Secondary Antibody [HRP (Horseradish Peroxidase)] |
| NB720-B | Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin] |
| H0000319-P01-10ug | Recombinant Human Apolipoprotein F Protein |
| H00003990-R01-20nmol | LIPC RNAi |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H0000319-B01P

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

