

Product Datasheet

TRIM23 Antibody (2G4) - Azide and BSA Free H00000373-M12

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00000373-M12

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00000373-M12



H00000373-M12

TRIM23 Antibody (2G4) - Azide and BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	2G4
Preservative	No Preservative
Isotype	IgG1 Kappa
Purity	Protein A purified
Buffer	In 1x PBS, pH 7.4


Product Description	
Description	Quality control test: Antibody Reactive Against Recombinant Protein.
Host	Mouse
Gene ID	373
Gene Symbol	TRIM23
Species	Human
Specificity/Sensitivity	TRIM23 (2G4)
Immunogen	TRIM23 (AAH22510, 1 a.a. ~ 574 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa. MATLVVNKLGAGVDSGRQGSRGTA VVKVLECGVCEDEVFSLQGDKVPRLLLCG HTVCHDCLTRLPLHGRAIRCPFDRQVTDLGDSGVWGLKKNFALLELLERLQNG PIGQYGAAEESIGISGESIIRCEDEAHLASVYCTVCATHLCSECSQVTHSTKTL AKHRRVPLADKPHEKTMCSQHQAIEFVCL EEGCQTSPLMCCVCKEYGKHQ GHKHSVLEPEANQIRASILDMAHCIRTFTEEISDYSRKLVGIVQHIEGGEQIVEDG IGMAHTEHVPGTAENARSCIRAYFYDLHETLCRQEEMALSVVDAHVREKLIWLR QQQEDMTILLSEVSAACLHCEKTLQQDDCRVVLAKQEITRLLET LQKQQQQFT EVADHIQLDASIPVTFTKDNRVHIGPKMEIRVVTLGLDGAGKTTILFKLKQDEFM QPIPTIGFNVETVEYKNLKFTIWDVGGKHKLRPLWKHYLNTQAVVFVVDSSHR DRISEAHSELAKLLTEKELRDALLIFANKQDVAGALSVEEITELLSLHKLCCGRS WYIQGCDARSGMGLYEGLDWLSRQLVAAGVLDVA
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details	
Applications	ELISA, Sandwich ELISA
Recommended Dilutions	ELISA, Sandwich ELISA
Application Notes	It has been used for ELISA.



Images

Sandwich ELISA: TRIM23 Antibody (2G4) [H00000373-M12] - Detection limit for recombinant GST tagged TRIM23 is 1 ng/ml as a capture antibody.



**Novus Biologicals USA**

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00000373-M12

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-43319-0.5mg	Mouse IgG1 Kappa Isotype Control (P3.6.2.8.1)
H00000373-P01-10ug	Recombinant Human TRIM23 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00000373-M12

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

