

# Product Datasheet

## Arylsulfatase A/ARSA Antibody

### H00000410-D01P

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/H00000410-D01P](http://www.novusbio.com/H00000410-D01P)

Updated 6/27/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/H00000410-D01P](http://www.novusbio.com/reviews/destination/H00000410-D01P)



**H00000410-D01P****Arylsulfatase A/ARSA Antibody**

| <b>Product Information</b> |                                                                                                                  |
|----------------------------|------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>           | 0.1 mg                                                                                                           |
| <b>Concentration</b>       | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| <b>Storage</b>             | Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.                                                     |
| <b>Clonality</b>           | Polyclonal                                                                                                       |
| <b>Preservative</b>        | No Preservative                                                                                                  |
| <b>Isotype</b>             | IgG                                                                                                              |
| <b>Purity</b>              | Protein A purified                                                                                               |
| <b>Buffer</b>              | PBS (pH 7.4)                                                                                                     |


| <b>Product Description</b>     |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
|--------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>             | Quality control test: Antibody reactive against mammalian transfected lysate.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Host</b>                    | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Gene ID</b>                 | 410                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Gene Symbol</b>             | ARSA                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Species</b>                 | Human                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Specificity/Sensitivity</b> | Reacts with arylsulfatase A.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Immunogen</b>               | ARSA (AAH14210.2, 1 a.a. - 507 a.a.) full-length human protein.<br>MGAPRSLLLALAAGLAVARPPNIVLIFADDLGYGDLGCYGHPSSTTPNLDQLAA<br>GGLRFTDFYVPVSLCTPSRAALLTGRLPVRMGMYPGVLPSSRGGLPLEEVT<br>AEVLAARGYLTMAGKWHLGVGPEGAFLPPHQGFHRFLGIPYSHDQGPCQNL<br>TCFPPATPCDGGCDQGLVPIPLLANSVEAQPPWLPGLEARYMAFAHDLMA<br>DADAQRQDRPFLLYYASHHHTYPQFSGQSFAERSGRGPFGLDLMELDA<br>AVGTLMTAIGDLGLLEETLVIFTADNGPETMRMSRGGCSGLLRGCGKGT<br>TYEGGVREPALAFWPGHIAPGVTHELASSLDLLPTLAALAGAPLPNV<br>TLDGFDLSPLLLGTGKSPRQSLFFYPSYPDEV<br>RGVFAVRSGKYKAHFFTQGS<br>AHSDDTADPACHASSSLTAHE<br>PPLLYDLSKDPGENYNLLGGVAGATPEVLQALKQLQLLKAQLDA<br>AVTFGPSQVARGEDPALQICCHPGCTPRPACCHCPDPA |
| <b>Notes</b>                   | This product is produced by and distributed for Abnova, a company based in Taiwan.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |

| <b>Product Application Details</b> |                                                                                                  |
|------------------------------------|--------------------------------------------------------------------------------------------------|
| <b>Applications</b>                | Western Blot                                                                                     |
| <b>Recommended Dilutions</b>       | Western Blot                                                                                     |
| <b>Application Notes</b>           | This antibody is reactive against transfected lysate in WB and as a detection antibody in ELISA. |



## Images

Western Blot: Arylsulfatase A/ARSA Antibody [H00000410-D01P] - Analysis of ARSA expression in transfected 293T cell line by ARSA polyclonal antibody. Lane 1: ARSA transfected lysate(53.60 KDa). Lane 2: Non-transfected lysate.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
[novus@novusbio.com](mailto:novus@novusbio.com)

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
[canada.inquires@bio-techne.com](mailto:canada.inquires@bio-techne.com)

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
[info.EMEA@bio-techne.com](mailto:info.EMEA@bio-techne.com)

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [technical@novusbio.com](mailto:technical@novusbio.com)  
Orders: [orders@novusbio.com](mailto:orders@novusbio.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis.  
Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/H00000410-D01P](http://www.novusbio.com/reviews/submit/H00000410-D01P)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

