

# Product Datasheet

## Bone marrow stromal cell antigen 1/CD157 Antibody Pair H00000683-AP11

Unit Size: 1 Set

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/H00000683-AP11](http://www.novusbio.com/H00000683-AP11)

Updated 6/27/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/H00000683-AP11](http://www.novusbio.com/reviews/destination/H00000683-AP11)



**H00000683-AP11****Bone marrow stromal cell antigen 1/CD157 Antibody Pair**

| <b>Product Information</b> |                                                                                                                       |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>           | 1 Set                                                                                                                 |
| <b>Concentration</b>       | Concentration of individual antibodies may be found on the vial label. If unlisted please contact technical services. |
| <b>Storage</b>             | Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.                                                          |
| <b>Clonality</b>           | Monoclonal                                                                                                            |
| <b>Preservative</b>        | No Preservative                                                                                                       |
| <b>Buffer</b>              | 1x PBS, pH 7.4                                                                                                        |

| <b>Product Description</b>  |                                                                                                                                                                                                                                                                                                                                                                                 |
|-----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>          | <p>Quality control test: Standard curve using recombinant protein (H00000683-P01) as an analyte. This antibody pair set comes with matched antibody pair to detect and quantify protein level of human BST1.</p> <p>Capture antibody: rabbit MaxPab affinity purified polyclonal anti-BST1 (100 ug)<br/> Detection antibody: mouse monoclonal anti-BST1 IgG2a Kappa (20 ug)</p> |
| <b>Gene ID</b>              | 683                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Gene Symbol</b>          | BST1                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Species</b>              | Human                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Notes</b>                | This product is produced by and distributed for Abnova, a company based in Taiwan.                                                                                                                                                                                                                                                                                              |
| <b>Standard Curve Range</b> | 3 ng/ml to 100 ng/ml.                                                                                                                                                                                                                                                                                                                                                           |
| <b>Assay Type</b>           | Sandwich ELISA                                                                                                                                                                                                                                                                                                                                                                  |

| <b>Product Application Details</b> |                                    |
|------------------------------------|------------------------------------|
| <b>Applications</b>                | Sandwich ELISA                     |
| <b>Recommended Dilutions</b>       | Sandwich ELISA 3 ng/ml - 100 ng/ml |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA

Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada

Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Antibody Pairs are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/H00000683-AP11](http://www.novusbio.com/reviews/submit/H00000683-AP11)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

