

Product Datasheet

AP1S1 Antibody - Azide and BSA Free H00001174-B02P

Unit Size: 0.05 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00001174-B02P

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00001174-B02P



H00001174-B02P

AP1S1 Antibody - Azide and BSA Free

Product Information

Unit Size	0.05 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	PBS (pH 7.4)

Product Description


Description	Quality control test: Antibody reactive against mammalian transfected lysate.
Host	Mouse
Gene ID	1174
Gene Symbol	AP1S1
Species	Human
Specificity/Sensitivity	adaptor-related protein complex 1, sigma 1 subunit
Immunogen	AP1S1 (NP_001274.1, 1 a.a. - 158 a.a.) full-length human protein. MMRFMLLFSRQGKLRQLQKWYLATSDKERKKMVRELMQVVLARKPKMCSFLE WRDLKVYKRYASLYFCCAIEGQDNELITLEIHRYVELLDKYFGSVCELDIIFNF EKAYFILDEFLMGGDVQDTSKKSVLKAIQADLLQEEDESPRSVLEEMGLA
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details

Applications	Western Blot
Recommended Dilutions	Western Blot
Application Notes	This antibody is reactive against mammalian transfected lysate in WB. It has been used for ELISA.

Images

Western Blot: AP1S1 Antibody [H00001174-B02P] - Analysis of AP1S1 expression in transfected 293T cell line by AP1S1 polyclonal antibody. Lane 1: AP1S1 transfected lysate(18.70 KDa). Lane 2: Non-transfected lysate.



**Novus Biologicals USA**

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00001174-B02P

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97019-5mg	Mouse IgG Isotype Control
H00001174-P01-10ug	Recombinant Human AP1S1 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00001174-B02P

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

